

DNA copy-number measurement of genome replication dynamics by high-throughput sequencing: the sort-seq, sync-seq and MFA-seq family

Dzmitry G. Batrakou^{1,2}, Carolin A. Müller^{1,2}, Rosemary H. C. Wilson^{1,2} and Conrad A. Nieduszynski^{1*}

Genome replication follows a defined temporal programme that can change during cellular differentiation and disease onset. DNA replication results in an increase in DNA copy number that can be measured by high-throughput sequencing. Here we present a protocol to determine genome replication dynamics using DNA copy-number measurements. Cell populations can be obtained in three variants of the method. First, sort-seq reveals the average replication dynamics across S phase in an unperturbed cell population; FACS is used to isolate replicating and non-replicating subpopulations from asynchronous cells. Second, sync-seq measures absolute replication time at specific points during S phase using a synchronized cell population. Third, marker frequency analysis can be used to reveal the average replication dynamics using copy-number analysis in any proliferating asynchronous cell culture. These approaches have been used to reveal genome replication dynamics in prokaryotes, archaea and a wide range of eukaryotes, including yeasts and mammalian cells. We have found this approach straightforward to apply to other organisms and highlight example studies from across the three domains of life. Here we present a *Saccharomyces cerevisiae* version of the protocol that can be performed in 7–10 d. It requires basic molecular and cellular biology skills, as well as a basic understanding of Unix and R.

Introduction

Complete and faithful DNA replication is essential for cell division, and, therefore, the development and maintenance of any organism. Bi-directional DNA replication forks initiate from origins of replication following their activation (or firing) during the synthesis or S phase of the cell cycle¹. In eukaryotes and archaea, efficient replication of the whole genome is ensured by the presence of multiple origins throughout the genome, not all of which will fire in each S phase². In organisms with multiple origins, the activity of individual origins or clusters of origins gives rise to a characteristic temporal pattern for genome replication³.

Fundamentally, there are three different types of methods used to study genome replication timing: DNA copy-number-based methods^{4,5}, which we describe here; incorporation of dense isotopes or nucleotide analogues that allow isolation of replicated DNA by density centrifugation or immunoprecipitation^{6,7}; and immunoprecipitation of replication fork proteins⁸. We compare the advantages and limitations of these methods below.

DNA copy number as a proxy for DNA replication

During S phase, the relative DNA copy number increases for each region as it is replicated; therefore, measuring DNA copy number can be used as a proxy for replication time^{4,9}. For eukaryotes, DNA copy number increases from 1N to 2N in haploids and 2N to 4N in diploids (with a few notable exceptions¹⁰; for example, a greater upper limit for organisms undergoing multiple, concurrent rounds of replication) as the cell passes through S phase. Regions of DNA that are copied very early during S phase will be at their maximum copy number for most of S phase, giving an average S phase copy number near to this maximum (e.g., 2N in haploids and 4N in diploids). Conversely, regions of

¹Sir William Dunn School of Pathology, University of Oxford, Oxford, UK. ²These authors contributed equally: Dzmitry G. Batrakou, Carolin A. Müller, Rosemary H. C. Wilson. *e-mail: conrad.nieduszynski@path.ox.ac.uk

DNA that are replicated very late have an average S phase copy number close to 1N in haploids and 2N in diploids. Regions of DNA that are replicated in the middle of S phase have an average S phase copy number that lies on a continuous scale between 1N and 2N in haploids or 2N and 4 N in diploids. Measuring the DNA copy number genome-wide at time points during S phase or in a population of cells that represent all stages of S phase can, therefore, serve as a proxy for the time at which a particular region was replicated. We note that, as these are population-based methods, these data are determined by a combination of both the efficiency of a replication origin being active (firing) within the population, or how likely an origin is to fire, and the time at which the origin fired within each cell. These two parameters can be interrogated mathematically⁹.

The protocol we describe here measures DNA copy number as a proxy for replication timing by high-throughput sequencing (HTS). Different approaches for obtaining samples give rise to a family of closely related methods: sort-seq, sync-seq and marker frequency analysis (MFA-seq). Our protocol uses *Saccharomyces cerevisiae* as an example, and we also show data from *Haloferax volcanii*. In addition, this family of methods has been used to determine replication timing for a wide variety of organisms (see Box 1). Each method compares DNA copy number between replicating and non-replicating samples, with the key difference being the manner in which the samples are obtained. In sort-seq, FACS is used to obtain replicating cells from S phase and a non-replicating sample^{11,12}. In sync-seq, cells are cell-cycle synchronized, and replicating samples are taken during S phase, with non-replicating samples taken before (or after) S phase⁹. MFA-seq measures DNA copy number in any rapidly proliferating cell sample normalized to a stationary phase (non-replicating) control^{13,14}. We detail consideration for which variant to use in the 'Choice of variant' section below.

Applications of the method

An increasing number of studies, from a range of laboratories, have used this family of methods to study replication dynamics in bacteria^{13,14}, archaea^{14,15} and eukaryotes^{16–18}. For the most straightforward application, these techniques allow replication timing to be determined for wild-type and genome-engineered cells, including cells with synthetic chromosomes^{19,20}. Replication dynamics can also be measured following specific cell treatments (for example, hydroxyurea to interrogate the effect of replication stress)⁴. When combined with experimentally induced genetic mutations, these techniques can be used to assess the impact of genetic inactivation of specific origins in *cis* (for example, to determine the effect of delaying the replication of specific genes or genomic loci)^{14,21–23}. These techniques can also be used to monitor the effect of naturally occurring (or disease-related) polymorphisms on replication dynamics¹⁷. Alternatively, the *trans* effect of genetic alterations can be investigated (for example, the genome-wide impact of altered replication protein expression either directly or relative to the number of active origins)^{22,24–27}. At a broader level, a number of studies have used these methods to examine the evolution of replication origin location and activity, giving further insights into the control of replication timing^{12,16,21}. The use of a standardized protocol that can be applied to a range of organisms is particularly important for these comparative genomic studies.

These techniques are also suitable for investigating non-typical DNA replication. For example, with slight adaptations, over-replication and under-replication can be investigated^{28,29}, as well as the identification of regions that replicate after S phase (e.g., during G2 phase)²⁴. In addition, these techniques have the potential to measure the DNA replication that occurs when cells repair DNA double-strand breaks by Break-Induced-Replication³⁰. Finally, data generated by this family of methods can be interrogated to reveal altered karyotypes, including large structural rearrangements^{13,14} and aneuploidy (see Step 61A(iii) and Box 2). Importantly, the protocol provides genome-wide datasets, which allow comparison with other datasets including protein binding sites (e.g., replication factors, transcription factors and chromatin marks), genetic elements (e.g., tRNA genes and autonomously replicating sequences) or structural features (e.g. G4 quadruplexes and topologically associating domains).

Alternative methods of measuring genome replication dynamics

A major advantage of the copy-number-based approaches described here is that they do not require genetic manipulations or special treatments, allowing analysis of unperturbed wild-type cells. However, the dynamic range is limited (e.g., from one to two for most eukaryotes), thus requiring an accurate quantification of copy number. This requires high coverage sequencing and can prove challenging for large genomes if high spatial resolution is required¹⁸. This requirement can be

Box 1 | Adapting for other organisms

Although the protocol we describe here uses *Saccharomyces cerevisiae* as the example organism, this family of methods can be applied to virtually any organism. Some adaptation to the protocol may be necessary for other organisms; nevertheless, we have found it straightforward to apply this family of methods across the three domains of life. Below, we provide general guidance for how to adapt the protocol and a list of organisms in which this protocol has been successfully applied.

Steps 1 and 2 (obtaining replicating and non-replicating cell populations) require organism-specific adaptation. As detailed above in Choice of variant, this is most straightforward for sort-seq and MFA-seq, where culture conditions should be optimized to maximize cell proliferation and thus obtain a high fraction of replicating cells. Note that for sort-seq, the most appropriate DNA stain and cell-sorting condition may need to be optimized; for example, we have found propidium iodide to be optimal when sorting cultured mammalian cells¹⁸. For sort-seq, the non-replicating samples are obtained by cell sorting of either the G1 or G2/M population, whichever is most numerous for the organism and growth conditions employed. For example, for mammalian cells, we have used sorted G1 cells as the non-replicating control sample. For sync-seq, a critical requirement is to have a cell population pass synchronously through S phase. While this generally requires an organism-specific protocol, elutriation has been successfully applied to a wide range of yeast species for which mating pheromone-based arrest and release approaches are not available^{16,48}.

The DNA extraction steps (3–17) should also be optimized to suit the organism of choice; we recommend following published protocols. However, we have observed an extraction bias in DNA samples that were prepared using DNA purification columns; therefore, we recommend phenol:chloroform extraction and ethanol precipitation. Once DNA has been obtained from the samples, library making and sequencing can be performed as described here (Steps 18–56), with consideration of the number of reads required based upon genome size (see Multiplex sequencing and sequence depth requirements in Experimental design). Finally, some organisms may benefit from optimization of the computational protocol, particularly, the read mapping steps (59 and 60). For example, we have found STAR⁴⁹ to be a more efficient read mapper than BOWTIE2 when working with the human genome¹⁸. More details on adapting the protocol to other organisms are available from the studies listed below.

Organism	Example reference
Fission yeast species	Sort-seq: <i>Schizosaccharomyces pombe</i> ⁴⁸ Sync-seq: G2 synchronization by centrifugal elutriation <i>Schizosaccharomyces pombe</i> ⁴⁸ G2 synchronization by centrifugal elutriation <i>Schizosaccharomyces pombe</i> , <i>Schizosaccharomyces octosporus</i> , <i>Schizosaccharomyces japonicus</i> ⁵⁰ Cdc25-22 temperature G2 arrest <i>Schizosaccharomyces pombe</i> ⁵¹
Budding yeast species	MFA-seq: <i>Lachanea fantastica</i> , <i>Lachanea meyersii</i> , <i>Lachanea dasiensis</i> , <i>Lachanea nothofagi</i> , <i>Lachanea waltii</i> , <i>Lachanea thermotolerans</i> , <i>Lachanea mirantina</i> , <i>Lachanea fermentati</i> , <i>Lachanea cidri</i> , <i>Lachanea kluyveri</i> ¹⁶ <i>Lachanea kluyveri</i> ⁵² <i>Lachanea waltii</i> ⁵³ Sort-seq: <i>Kluyveromyces lactis</i> , <i>Lachanea kluyveri</i> , <i>Candida glabrata</i> , <i>Naumovozyma castellii</i> , <i>Tetrapispora blattae</i> , <i>Zygosaccharomyces rouxii</i> ²¹ <i>Saccharomyces paradoxus</i> , <i>Saccharomyces arboricolus</i> , <i>Saccharomyces bayanus</i> ¹² Sync-seq: G1 synchronization by centrifugal elutriation <i>Lachanea fantastica</i> , <i>Lachanea meyersii</i> , <i>Lachanea dasiensis</i> , <i>Lachanea nothofagi</i> , <i>Lachanea waltii</i> , <i>Lachanea thermotolerans</i> , <i>Lachanea mirantina</i> , <i>Lachanea fermentati</i> , <i>Lachanea cidri</i> , <i>Lachanea kluyveri</i> ¹⁶
Other unicellular eukaryotes	Sort-seq: <i>Trypanosoma brucei</i> ^{41,54,55} <i>Leishmania major</i> , <i>Leishmania Mexicana</i> ⁵⁶
Archaea	MFA-seq: <i>Sulfolobus acidocaldarius</i> , <i>Sulfolobus solfataricus</i> ⁵⁷ <i>Haloferax volcanii</i> ¹⁴ Sync-seq: Acetic acid arrest <i>Sulfolobus acidocaldarius</i> ⁵⁷
Bacteria	MFA-seq: <i>Escherichia coli</i> ²⁶
<i>Homo sapiens</i>	MFA-seq: 161 individuals sequenced by the 1000 Genomes Project ¹⁷ Sort-seq: CO202 (microarray) ⁵⁸ HeLa ¹⁸ Six lymphoblastoid cell lines ¹⁷

Box 2 | Genotypic information from replication profiles

Genome-wide DNA copy-number measurements can reveal (unanticipated) genotypic variations.

▲ CRITICAL Visualize the non-replicating sample without normalization.

Aneuploidy: Most cells contain complete chromosome sets, with each chromosome present at the same copy number. Ploidy differences result in individual chromosomes with higher (or lower) copy number than the remaining genome (see Fig. 6a,c). Stable aneuploid karyotypes result in integer copy numbers, whereas fractional increases (or decreases) indicate heterogeneity in the population.

Segmental duplications: Chromosomal regions with a copy number greater than the flanking sequences and other chromosomes indicate localized sequence amplifications (see Fig. 6b,d).

Chromosomal rearrangements: Adjacent chromosomal sequences replicate at very similar times in S phase, giving rise to smooth replication profiles. Steep jumps in the copy number of adjacent sequences often indicate chromosomal rearrangements relative to the reference sequence¹⁴. The karyotype of the cell population should be analyzed further.

mitigated by alternative approaches that utilize labelling of nascent DNA with nucleotide analogues, followed by enrichment of newly synthesized DNA. A variant of nascent DNA labelling and separation is the dense isotope transfer method similar to the Meselson–Stahl experiment; however, it is expensive and labor-intensive and has separation biases³¹. More recent variants include 5-bromo-2'-deoxyuridine-immunoprecipitation-sequencing (BrdU-IP-seq)⁷, 5-ethynyl-2'-deoxyuridine-sequencing in hydroxyurea (EdUseq-HU)^{32,33} and repli-seq³⁴. In these methods, nascent DNA containing an incorporated nucleotide analogue, such as BrdU or EdU, is enriched through immunoprecipitation or click chemistry and identified by microarrays or HTS. The enrichment for nascent DNA increases the dynamic range and thus makes these approaches particularly useful for organisms with large genomes. However, this qualitative separation of nascent from parental DNA comes at the expense of quantitative information, particularly biases from the enrichment, and, consequently, the measurements are expressed in arbitrary units.

An alternative method for studying replication timing is through ChIP-seq of replication fork proteins (e.g., go-ichi-ni-san (GINS) components)⁸. However, this approach is more technically challenging, requiring cell synchronization and either tagging of proteins by genetic alteration or the availability of a suitable antibody to the endogenous protein.

All of these approaches to measure genome replication dynamics can be complemented with methods that identify replication origins^{35,36}, determine replication fork direction^{37–39} or investigate the response to replication stress⁴⁰.

Limitations

The family of methods described here are powerful techniques to study replication timing in various organisms. Nevertheless, there are, as always, some limitations of the techniques. Firstly, mapping of the high-throughput sequence data requires a reference genome, which is not available for some organisms. Second, the sort-seq and MFA-seq variants produce replication profiles by measuring the relative copy number of all genomic loci; i.e., there is no absolute measure of replication time. This is in contrast to sync-seq, where DNA copy number is measured at specific time points to allow each locus to be assigned a median replication time relative to a cell cycle landmark, such as the start of DNA synthesis. Third, as mentioned above, the small dynamic range between one and two in eukaryotes requires an accurate quantification of DNA copy number, which can limit spatial resolution in organisms with large genomes. In this case, if only a few regions of the genome are of interest, loci-specific copy number quantitative polymerase chain reaction (qPCR)⁴¹ or droplet digital polymerase chain reaction (ddPCR)¹⁸ methods can be used. Such approaches can also be used to test for allele-specific replication times^{17,18}. Fourth, a non-replicating population provides a control to minimize sequence biases. However, replication profiles can be generated without a non-replicating control sample (Step 61A(v)), but they may be noisier. Finally, it is also important to reiterate that this protocol generates data from a population of cells and as such reports the population-average replication time. If relevant for the user, copy number approaches have recently been adapted to single cells⁴², and labelled nascent DNA approaches have recently been adapted to single molecules⁴³.

Overview of the procedure

The first section of the protocol (Steps 1 and 2, options A, B or C) describes the experimental setup and cell growth for each of the three variants of the family: sort-seq, sync-seq and MFA-seq




Fig. 1 | Workflow overview. For each variant, the protocol starts with obtaining replicating (cyan) and non-replicating (gray) cells (Steps 1 and 2). Sort-seq uses FACS to enrich cells in S phase or late G2 phase as indicated. Sync-seq uses cell cycle arrest and release to collect a non-replicating G1 sample and replicating samples at time points during S phase. In MFA-seq, log and stationary phase cultures provide the replicating and non-replicating samples, respectively. Criteria for selecting the most appropriate variant are discussed in Choice of variant. The DNA isolated from the samples (Steps 3–17) is prepared for HTS using short-read technology (Steps 18–56). The resulting reads are aligned and summed in genomic bins (Steps 57–60); example replicating and non-replicating data are shown with the whole genome median read count represented by a horizontal pink line. Finally, the replicating data can be normalized using the non-replicating data to generate replication timing profiles (Step 61); an example replication profile is shown with smoothing in orange.

(see Fig. 1). The variants converge for the following experimental sections of the protocol: DNA extraction (Steps 3–17), library preparation (Steps 18–52) and HTS (Steps 53–56). The last section of the protocol details data analysis, firstly aligning and processing of sequencing data (Steps 57–60) and then normalization and visualization using Repliscope, our R package available from the Comprehensive R Archive Network (CRAN), either using the command line (Step 61, option A) or a built-in graphical user interface (Step 61, option B). We cover validation and quality-control steps in both the experimental and data analysis parts of the protocol. At the end of the protocol, users will have generated replication timing profiles for their strain(s) or conditions of interest and compared them as appropriate. Finally, we suggest various bespoke analyses that users can consider undertaking with their data if desired.

Experimental design

Choice of variant

The choice of variant depends on both the organism and experimental question. Sync-seq provides the highest temporal and spatial resolution of the three variants; it is the only variant that provides an absolute measure of replication time. However, this requires S phase cell cycle synchronization from

Box 3 | Normalization

All three variants of the method rely on a non-replicating control sample to normalize for coverage biases from the genomic sequence context, including non-unique sequences that are discarded during the analysis.

Sort-seq takes advantage of a sorted non-replicating fraction, either G1 or G2/M.

Sync-seq uses an arrested, unreleased G1 sample.

MFA-seq uses a stationary phase culture as a non-replicating sample.

While the above ways of obtaining non-replicative samples are commonly used with the respective methods, there is no strict requirement for a particular one. For example, an arrested sample may be used as a non-replicative control in either MFA-seq or sort-seq.

▲ CRITICAL Visualize the non-replicating sample without normalization to make sure that the read counts profile is even across all regions.

The first normalization step is to calculate ratios between read counts of a replicating sample and a non-replicating sample, adjusted for the differences in total read number between the two samples. The resulting ratios have a mean of 1.

The second normalization involves multiplying the ratio values by a normalization factor to fit the data to a biologically relevant relative copy-number scale (in most cases from 1 to 2). This is calculated differently for sync-seq, sort-seq and MFA-seq.

Sync-seq samples are normalized to their bulk genome content measured by flow cytometry. For accuracy, the median DNA content for every time point should be plotted and fitted to a sigmoid function, before determining the percentage of genome replication for relevant time points. This can be done for experimental flow cytometry data using our online tool: <https://dnareplab.shinyapps.io/fitsigmoid/>. The normalization factor for individual time points is then calculated as $1 + [\text{fraction replicated}]$. For example, for a sample with 50% bulk genome replicated, the normalization factor would be 1.5.

For sort-seq samples, relative copy-number values range between 1 and 2. Therefore, the normalization factor is calculated to maximize the number of data points lying within this range. Normalization factors of 1.3–1.5 are usually appropriate. This is implemented in the ‘Auto normalisation’ function of the Repliscope package.

For samples containing both replicating and non-replicating cell populations (i.e., MFA-seq), a similar approach can be used, but with a different maximum value for relative copy number. The value can be calculated as $1 + [\text{proportion of replicating cells}]$. For example, an asynchronous culture of *Saccharomyces cerevisiae* contains ~25% replicating cells. Thus, the expected relative copy number range is 1 to 1.25.

which samples can be taken as genome replication proceeds. Sort-seq alleviates the requirement for cell cycle synchronization but produces relative rather than absolute replication timing data. The main requirement for sort-seq is that replicating cells are clearly identifiable in a DNA-content flow cytometry profile, since this allows a representative replicating cell population to be obtained by FACS. Cell sorting may require optimization of DNA staining conditions and cell sorter setup for the organism of interest; this may be problematic for organisms with a variable ploidy or that undertake cell septation during S phase. MFA-seq is technically the most straightforward method to apply to any organism, as it solely requires proliferating and stationary phase samples. However, the resulting dynamic range is dependent on the fraction of replicating cells within the population. Therefore, MFA-seq is particularly well suited to organisms for which the other variants are not technically applicable (e.g., those that undergo concurrent rounds of replication, cannot be synchronized or have unstable ploidy)¹⁴. For organisms where multiple variants could be used, MFA-seq typically generates the smallest dynamic range and thus will be of lower resolution. Box 1 provides further information for adapting this protocol to other organisms, as well as a list of published uses of this family of methods across the three domains of life.

Controls

Replicating samples should be normalized (in almost all cases) to a non-replicating sample to control for sequencing bias across the genome. The key requirement for the non-replicating control is that the whole genome is at a uniform copy number. Typically, the non-replicating control is DNA from G1 or G2 phase cells for sort-seq, DNA from cells arrested in G1 phase for sync-seq and DNA from stationary cultures for MFA-seq. Sometimes the use of alternative non-replicating controls may be required; we detail these and other considerations pertaining to the users’ choice of normalization control in Box 3. The non-replicating control sample should be prepared in parallel and sequenced on the same run as replicating samples.

Multiplex sequencing and sequence depth requirements

By using multiplexing barcodes in Step 42 of the protocol, users can combine multiple samples on the same sequencing run. Reads are then demultiplexed at the end of the run using the differing barcode sequences. However, before multiplexing samples, users must consider the optimal number of reads per sample, which depends heavily on the genome size of the organism of interest and the chosen bin

size. A simple formula to calculate the number of multiplexed samples per run is $N = (t \times b)/(g \times d)$, where t is the estimated total yield of the run (reads), b is bin size (bp), g is haploid genome size (bp) and d is desired sequencing depth per bin (reads). We find that a depth of 1,000–2,000 reads/bin gives cost-efficient sequencing depth without compromising the precision of the data. By default, we perform single-end sequencing, although the computational pipeline described here can be applied to paired-end sequencing data. A NextSeq run typically yields ≥ 400 million single-end reads. Therefore, for samples from *Saccharomyces cerevisiae* (~13-Mbp genome), we combine up to 24 samples per NextSeq run for a sequencing depth of $>1,250$ reads/1-kb bin. In contrast, for a human cell line, we recently combined two samples on a NextSeq run to obtain an average depth of 2,095 reads/50-kb bin¹⁸. For more discussion of optimal read number for copy number-based methods, see ref. ⁴.

Required expertise

The protocol described here requires standard microbiology and molecular biology techniques. In addition, for in-house HTS, the experimentalist should be able to operate the sequencer (e.g., an Illumina NextSeq machine); an alternative is to use an HTS facility or company. The sort-seq protocol requires use of a FACS machine, either directly or through a facility. The data analyses require a basic familiarity with the Linux command line and the R programming language. For those with limited R experience, a web application provides a graphical user interface to R that allows the user to perform the standard data analyses described in this protocol. For more experienced users, the R Repliscope package permits an equivalent analysis using a command line interface and provides a starting point for more sophisticated custom analyses.

Materials

Reagents

Biological materials

Yeast strain of interest, taken from a fresh YPAD plate or appropriate Synthetic Complete (SC) medium (specific to *Saccharomyces cerevisiae*). For the example data, we used *Saccharomyces cerevisiae* yeast strains of a W303 genetic background (e.g., strain BMA64-1A from Euroscarf) (Figs. 1–6). Figure 5 also shows example MFA data for *H. volcanii*¹⁴ ▲ **CRITICAL** Yeast strains of mating type *a* are required for sync-seq ▲ **CRITICAL** For sync-seq, we used a *BAR1* strain. *BAR1* encodes a secreted protease that degrades alpha factor; therefore, with *bar1* strains, less alpha factor is required to arrest cells.

Chemicals

- 2-mercaptoethanol (Sigma-Aldrich, cat. no. M6250) ▲ **CRITICAL** Toxic and irritant; avoid inhalation and wear personal protective equipment (PPE).
- Agencourt AMPure XP magnetic beads (Beckman Coulter, cat. no. A63881)
- Alpha factor (specific to *Saccharomyces cerevisiae*; GenScript, cat. no. RP01002)
- Chloroform (Sigma, cat. no. 24216) ▲ **CRITICAL** Harmful and irritant; avoid inhalation and wear PPE.
- D-Sorbitol (Sigma-Aldrich, cat. no. S1876)
- Ethanol, absolute (Sigma-Aldrich, cat. no. 32221) ▲ **CRITICAL** Flammable.
- EDTA (0.5 M sterile solution pH 8.0; VWR, cat. no. E177-500)
- HCl (Sigma-Aldrich, cat. no. 30721) ▲ **CRITICAL** Corrosive; avoid inhalation and wear PPE.
- NEBNext Library Quant Kit for Illumina (NEB, cat. no. E7630S)
- NEBNext Multiplex Oligos for Illumina (Index Primers Set 1) (NEB, cat. no. E7335S)
- NEBNext Multiplex Oligos for Illumina (Index Primers Set 2) (NEB, cat. no. E7500S)
- NEBNext Ultra II DNA Library Prep Kit for Illumina (NEB, cat. no. E7645S)
- NextSeq 500/550 High Output Kit v2.5 (75 Cycles) (Illumina, cat. no. 20024906)
- Phenol:chloroform:isoamyl (Sigma, cat. no. 77617) ▲ **CRITICAL** Phenol is toxic and corrosive. Chloroform is harmful and an irritant. Wear PPE and avoid inhalation.
- Pronase (specific to *Saccharomyces cerevisiae*; Calbiochem/Merck, cat. no. 53702-250KU)
- Proteinase K (MP Biomedicals, cat. no. 193504) ▲ **CRITICAL** We have found that proteinase K from other suppliers can produce artefacts in flow cytometry DNA content analysis.
- RNase A (MP Biomedicals, cat. no. 219398050)
- Sodium dodecyl sulfate, 20% (wt/vol) solution (Fisher Scientific UK; cat. no. 10607633)
- Sodium azide (Sigma-Aldrich, cat. no. S8032) ▲ **CRITICAL** Toxic; handle with care and wear PPE.
- Sodium chloride (Sigma-Aldrich, cat. no. S9888)

- SYTOX Green Nucleic Acid Stain (5 mM solution; Life Technologies, cat. no. S7020)
- ▲ **CRITICAL** SYTOX Green generates superior flow cytometry DNA content profiles compared to propidium iodide for *Saccharomyces cerevisiae*⁴⁴.
- Tris base (Trizma base; Sigma-Aldrich, cat. no. T6066)
- Tri-sodium citrate (Fisher, cat. no. S/3320/53) ▲ **CRITICAL** We have found that tri-sodium citrate from other suppliers can produce artefacts in flow cytometry DNA content analysis.
- Qubit double-stranded DNA (dsDNA) High-Sensitivity (HS) Assay Kit (Life Technologies, cat. no. Q32854)
- TapeStation High Sensitivity D1000 Reagents (Agilent, cat. no. 5067-5585)
- TapeStation High Sensitivity D1000 Screentape (Agilent, cat. no. 5067-5584)
- YPAD broth or appropriate SC broth (specific to *Saccharomyces cerevisiae*; ForMedium, cat. no. CCM1005)
- YPAD agar or appropriate SC agar (specific to *Saccharomyces cerevisiae*; ForMedium, cat. no. CCM0510)
- Zymolyase T20 (specific to *Saccharomyces cerevisiae*; *Arthrobacter luteus*; Seikagaku, cat. no. 120491)

Equipment

- 1l glass, conical (e.g., VWR, cat. no. 214-1134)
- 50-ml centrifuge tubes (e.g., Greiner Bio One, cat. no. 227261)
- 100-ml glass, conical (e.g., Fisher, cat. no. FB33131)
- 4200 TapeStation System with personal computer (PC; Agilent, cat. no. G2991AA) or alternate DNA fragment length measurement (e.g., PAGE gels)
- Autoclave (e.g., Rodwell, Phoenix)
- Automatic pipette controller (e.g., Star Lab, ErgoOne Fast)
- Balance (e.g., Sartorius, cat. no. BP1200)
- Basespace (Illumina): set up a Basespace account (<https://basespace.illumina.com/home/index>)
- Bioruptor Standard (Diagenode, cat. no. UCD-200) or alternate sonicator suitable for DNA shearing
- Benchtop centrifuge for 15-ml and 50-ml centrifuge tubes (e.g., Thermo Scientific, Heraeus Multifuge X3R)
- Cuvettes (e.g., Fisher Scientific, Semi-micro cuvettes, 2.5 ml, cat. no. 10781791)
- Filter (e.g., Sartorius, cat. no. 16555-K)
- Freezer, −20 °C
- Refrigerator, 4 °C
- Fume hood
- Hemocytometer (e.g., Hawksley, cat. no. Z4+AC1000)
- Heat block (e.g., Techne, Dri-block DB-2D)
- Ice bucket
- Ice machine (e.g., Scotsman, cat. no. AF206)
- Incubator shaker (e.g., New Brunswick Scientific, Innova 40)
- Light microscope (e.g., Olympus, cat. no. CKX41)
- Magnetic microcentrifuge tube rack (e.g., Invitrogen, DynaMag-2, cat. no. 12321D)
- Microcentrifuge capable of chilling to 4 °C (e.g., Eppendorf, Centrifuge 5424R)
- Microcentrifuge tubes (e.g., Star Lab, cat. no. E1415-1500)
- Micropipettes (e.g., Star Lab, ErgoOne series)
- Micropipette tips (e.g., Star Lab, TipOne series)
- Microscope slides (e.g., Thermo Scientific, Superfrost Plus cat. no. J1800AMNZ)
- Milli-Q water purifier (e.g., Veolia, Purelab Flex)
- NanoDrop spectrophotometer (NanoDrop, cat. no. ND-1000) or another DNA concentration quantification
- NextSeq 500 (Illumina)
- PCR tubes or strips (e.g., Greiner Bio One, cat. no. 683201)
- Petri dishes (e.g., Sarstedt, cat. no. 82.1473)
- pH meter (e.g., Jenway, cat. no. 3510)
- qPCR tubes (e.g., StarLab, Rotor-Gene Style 4-Strip Tubes and Caps, 0.1 ml, cat. no. I1402-0400)
- Quantitative PCR machine (e.g., Qiagen, Rotor-Gene Q 2plex Platform, cat. no. 9001550)
- Qubit Assay Tubes (Invitrogen, cat. no. Q32856)
- Qubit Fluorometer (e.g., Invitrogen, cat. no. Q32857) or another high sensitivity DNA concentration quantification

- Rotating wheel (e.g., Stuart, Rotator SB3)
- Screw-capped microcentrifuge tubes (e.g., Star Lab, cat. no. E1415-2231)
- Serological pipettes (e.g., Greiner Bio One)
- Shaking heat block for microcentrifuge tubes (e.g., Eppendorf, Thermomixer Comfort)
- Spectrophotometer (e.g., ThermoScientific, Biomate 3, cat. no. 335905P-02)
- Syringe (e.g., Terumo, cat. no. SS+20ES1)
- Thermal Cycler (e.g., Biorad, Tetrad 2)
- Tip Sonicator (e.g., Branson, Digital Sonifier)
- Timer (e.g., Fisher Scientific, Traceable Nano Timer)
- Vortex (e.g., Fisher Scientific, TopMix, cat. no. FB15013)
- Water bath shaker (e.g., New Brunswick Scientific, Innova 3100)

Flow cytometry

- BD LSRFortessa X-20 flow cytometer (Becton Dickinson) equipped with 488-nm laser (Coherent), 505-nm long-pass and 530/30-nm band-pass filters
- PC running FACSDiva software
- FlowJo 10.5.3 software (FlowJo) or alternative software to analyze flow cytometry data
- MoFlo XPD cell sorter (Beckman-Coulter) equipped with 488-nm laser (Sapphire, Coherent), 505-nm short-pass and 530/40-nm band-pass filters or another comparable cell sorter
- PC running Summit software

Software and hardware

- Workstation with 64-bit processor and ≥ 4 GB of RAM and latest Ubuntu operating system (OS) with the following installed software: FastQC (<https://www.bioinformatics.babraham.ac.uk/projects/fastqc/>), Bowtie 2⁴⁵ (<http://bowtie-bio.sourceforge.net/bowtie2/index.shtml>), SAMtools⁴⁶ (<http://www.htslib.org/>), bedtools⁴⁷ (<https://bedtools.readthedocs.io/en/latest/>), Picard (<https://broadinstitute.github.io/picard/>), R (<https://www.r-project.org/>) and R package Repliscope (<https://cran.r-project.org/web/packages/Repliscope/>)
- Additionally, a web browser is required for the interactive analysis using the app. We recommend Google Chrome (<https://www.google.com/chrome>) or the RStudio built-in browser (<https://www.rstudio.com/>)

The analysis presented here has been performed with the following software versions: Ubuntu 18.04.01 LTS, fastqc 0.11.5, bowtie2 2.3.4.1, samtools 1.7, bedtools 2.26.0, picard-tools 2.8.1, R 3.4.4, Repliscope 1.0.0, ggplot2 3.0.0, shiny 1.1.0, colourpicker 1.0 and Google Chrome 70.0.3538.77 (available as a bootable Ubuntu disc image at <https://ln1.path.ox.ac.uk/groups/nieduszynski/Replibuntu/Replibuntu-18.04.0-amd64.iso.gz>)

Reagent setup

▲ CRITICAL All solutions are prepared with Milli-Q water and stored at room temperature (RT, 20 °C) for ≤ 12 months unless otherwise indicated.

3 mM alpha factor

Dissolve aseptically 10 mg of alpha factor in 2 ml of sterile water. Store in 0.5-ml aliquots at -80 °C. Avoid excessive freeze/thawing.

1 M D-sorbitol

Dissolve 18.2 g of D-sorbitol in 100 ml of water. Sterilize by autoclaving.

2 M D-sorbitol

Dissolve 36.4 g of D-sorbitol in 100 ml of water. Sterilize by autoclaving.

0.2 M EDTA, 0.1% (wt/vol) sodium azide

Combine 40 ml of 0.5 M EDTA pH 8.0, 1 ml of 10% (wt/vol) sodium azide and 59 ml of water.

70% (vol/vol) ethanol

Dilute 70 ml of 100% (vol/vol) ethanol to 100 ml with deionized H₂O.

80% (vol/vol) ethanol

Dilute 80 ml of 100% (vol/vol) ethanol to 100 ml with deionized H₂O. **▲ CRITICAL** For AMPure XP bead cleanup steps, it is important that 80% (vol/vol) ethanol is made fresh on the day.

20 mg/ml pronase

Dissolve 200 mg of pronase in 10 ml of water. Filter-sterilize using a 0.2-µm filter and syringe and store in 1-ml aliquots at −20 °C.

20 mg/ml proteinase K

Dissolve the whole 100-mg bottle in 5 ml of water and store in 0.5-ml aliquots at −20 °C.

10 mg/ml RNase A

Dissolve the whole 50-mg bottle in 5 ml of water and store in 1-ml aliquots at −20 °C.

10% (wt/vol) sodium azide

Dissolve 10 g of sodium azide in 100 ml of water. **▲ CRITICAL** Sodium azide is toxic. Wear PPE.

3 M sodium chloride

Dissolve 17.53 g of sodium chloride in 100 ml of water.

5 M sodium chloride

Dissolve 29.2 g of sodium chloride in 100 ml of water.

1 M sodium citrate pH 7.4

Dissolve 294.1 g of sodium citrate in 500 ml of water. Adjust the pH to 7.4 with HCl. Add additional water to a total volume of 1 l. Sterilize by autoclaving.

50 mM sodium citrate pH 7.4

Dilute 1 M sodium citrate pH 7.4 to 50 mM (1:20) with water. Filter to remove aggregates immediately before use by using a 0.45-µm filter.

Solution 1

Component	Amount (ml)	Final concentration
Sorbitol (2 M)	30	1.2 M
Tris-HCl pH 8.0 (1 M)	10	0.2 M
EDTA pH 8.0 (0.5 M)	2	0.02 M
Water	To 50 ml total	-
Total	50	-

Solution 2

Component	Amount (ml)	Final concentration
Tris-HCl pH 8.0 (1 M)	2.5	0.05 M
EDTA pH 8.0 (0.5 M)	10	0.1 M
NaCl (3 M)	1.7	0.1 M
SDS (20% wt/vol)	1.25	0.5% (wt/vol)
Water	To 50 ml total	-
Total	50	-

1× SYTOX Green in 50 mM sodium citrate

Dilute SYTOX Green (1:5,000) in 50 mM sodium citrate pH 7.4.

1× SYTOX Green, 0.1% (wt/vol) sodium azide in 50 mM sodium citrate

Dilute SYTOX Green (1:5,000) and 10% (wt/vol) sodium azide (1:100) in 50 mM sodium citrate pH 7.4.

1 M Tris-HCl pH 8.0

Dissolve 60.6 g of Tris base in 400 ml of water. Adjust the pH to 8.0 with HCl. Add additional water to a final volume of 500 ml. Sterilize by autoclaving.

10 mM Tris-HCl pH 8.0

Dilute 1 M Tris-HCl pH 8.0 (1:100) in water.

10 mM Tris-HCl pH 8.0, 10 mM NaCl

Dilute 1 M Tris-HCl pH 8.0 (1:100) and 5 M NaCl (1:500) in water.

TE buffer (10 mM Tris, 1 mM EDTA)

Dilute 1 M Tris-HCl pH 8.0 (1:100) and 0.5 M EDTA pH 8.0 (1:500) in water.

YPAD agar

Dissolve 56 g of YPAD agar powder in 800 ml of deionized water. Sterilize by autoclaving and cool to 60 °C. Pour 20 ml aseptically into plates and let the agar set for 10–20 min. Store hermetically sealed at 4 °C.

YPAD broth

Dissolve 40 g of YPAD broth powder in 800 ml of deionized water. Sterilize by autoclaving.

10 mg/ml zymolyase T20

Dissolve 100 mg in 10 ml of water and store in 1-ml aliquots at –20 °C.

Equipment setup**Sequencing analysis software**

All the software can be installed from Ubuntu's Universe repository using the following commands:

```
$ sudo add-apt-repository universe
$ sudo apt update
$ sudo apt install fastqc bowtie2 samtools bedtools picard-tools r-base
```

Repliscope

Installing Repliscope from CRAN will also install its dependencies: ggplot2, shiny and colourpicker. Repliscope can be installed with the following commands:

```
$ sudo R
> install.packages('Repliscope')
```

Procedure**Growing primary cell culture ● Timing 17 h (hands-on time: 10 min)**

- 1 Inoculate a single colony of the yeast strain into 5 ml of liquid medium (YPAD or appropriate SC selection) in a 50-ml centrifuge tube and incubate overnight in a shaking incubator at 30 °C (or 23 °C for sync-seq) with 200 rpm shaking.

Obtaining replicating and non-replicating cells ● Timing 1.5–3 d (see Step 2A, B or C)

- 2 Obtain samples of replicating cells and a non-replicating control. Follow option A for sort-seq (FACS enrichment of cells in S or G2 phase), option B for sync-seq (replicating samples from a synchronous S phase population) or option C for MFA-seq (proliferating and stationary phase cells).

(A) Sort-seq ● Timing 2 d (hands-on time: 4 h)

- (i) Measure cell density by spectrophotometry or a hemocytometer and inoculate 3×10^7 cells from the overnight culture into 25 ml of YPAD medium in a 100-ml conical flask; allow to




Fig. 2 | Cell sorting gates. Cell sorting for S and G2 phases from *Saccharomyces cerevisiae*. Cells were labelled with SYTOX Green as described in Step 2A(viii). **a**, A typical ungated plot showing DNA content signal (FL1) Height/Area. Cells are gated to remove doublets and debris (Step 2A(x)). **b**, Histogram shows DNA content after gating. Sorting gates for S and G2 phases of the cell cycle, as well as the ‘ALL’ control fraction, are indicated (Step 2A(xi and xii)). **c**, Typical flow cytometry profiles of the DNA content from sorted cell populations to check sort purity (Step 2A(xviii)).

grow to $0.7\text{--}1.5 \times 10^7$ cells/ml at 30 °C in a shaking incubator with 200 rpm shaking (~4.5 h). Check cell density by spectrophotometry or with a hemocytometer. For a wild-type W303 strain, the required cell density corresponds to an optical density (measured at 660 nm) of 0.5–0.9. The conversion between optical density and cell density is non-linear and strain dependent; therefore, we advise use of a hemocytometer to determine this relationship.

- (ii) At a density of $0.7\text{--}1.5 \times 10^7$ cells/ml, transfer culture into a 50-ml centrifuge tube and harvest cells by centrifugation at 2,000g for 5 min at 20 °C.
- (iii) Discard the medium. Resuspend the cell pellet in 25 ml of Milli-Q water. Pellet cells by centrifugation at 2,000g for 5 min at RT and discard the supernatant. Repeat this wash and centrifuge again to pellet cells.
- (iv) Fix cells by resuspending in 10 ml of 70% (vol/vol) ethanol and incubate for ≥ 1 h at RT (20 °C) or store at 4 °C overnight.

■ **PAUSE POINT** Cells can be stored at 4 °C for ≤ 1 year.

- (v) Pellet cells by centrifugation at 2,000g for 5 min at RT and discard the supernatant. To wash the cell pellet with sodium citrate, resuspend cells in 25 ml of freshly filtered 50 mM sodium citrate and pellet by centrifugation at 2,000g for 5 min at RT. Repeat this wash, pellet by centrifugation as before and resuspend cells in 1 ml of 50 mM sodium citrate.
- (vi) Add 250 μ l of 10 mg/ml RNase A and incubate for 1 h at 37 °C in a water bath.
- (vii) Add 100 μ l of 20 mg/ml proteinase K and incubate for 1 h at 55 °C in a water bath. Gently tap the tube every 20 min to resuspend cells. Pellet cells by centrifugation at 2,000g for 5 min at RT.
- (viii) Stain DNA by adding 2.5 ml of 50 mM sodium citrate, 25 μ l of 10% (wt/vol) sodium azide and 5 μ l of SYTOX Green. Incubate in the dark for ≥ 1 h at RT or at 4 °C overnight.

■ **PAUSE POINT** SYTOX Green stained cells can be stored in the dark at 4 °C for ≤ 1 year.

- (ix) Immediately prior to FACS, visually inspect the cell aliquot by light microscopy. If cells are sticking to each other, sonicate the sample twice for 3 s using a tip sonicator set to 20 W of sonic power or vortex vigorously until cells are separate.

▲ **CRITICAL STEP** Avoid excessive sonication, as it damages cells and causes debris in the sample.

- (x) On a MoFlo XFP cell sorter (or another comparable cell sorter), use the signal height and area information from the 488-nm laser (as a measure of DNA content) to select the desired population of cells (S and G2 phase, or ‘ALL’ population). Exclude cell doublets or debris, as shown in Fig. 2a. Create a histogram with cell count on the y axis and signal area of DNA content on the x axis.

? TROUBLESHOOTING

- (xi) Set the sorting gate to include single cells from all cell cycle stages and collect 300,000 cells into a 50-ml centrifuge tube as the ‘ALL’ control sample, as shown in Fig. 2b.
- (xii) Adjust sorting gates to two populations representing S and late G2 phase cells as shown in Fig. 2b. For both fractions, collect ≥ 5 million cells into separate 50-ml centrifuge tubes.

- (xiii) Add 1/3 volume of 100% (vol/vol) ethanol to the FACS-enriched cells and mix by inversion.

■ **PAUSE POINT** Samples can be stored at 4 °C overnight.

- (xiv) For each sorted sample (i.e., there are now three sorted samples per strain: 'ALL', 'S' and 'G2'), the sorted cells are pelleted in a microcentrifuge tube. For each sorted sample, transfer 1.5 ml into a fresh 1.5-ml microcentrifuge tube and collect cells by centrifugation at 20,000g for 10 min at RT. Carefully remove and discard 1.4 ml of the supernatant without disturbing the cell pellet.
 - (xv) Transfer another 1.4 ml of the same sample into the same microcentrifuge tube, centrifuge and remove the supernatant as before. Repeat these steps until all FACS-enriched cells have been pelleted into one microcentrifuge tube per sorted sample. Carefully remove the supernatant after the final centrifugation without disturbing the cell pellet.
 - (xvi) Resuspend the S and G2 phase samples in 500 µl of Milli-Q water and transfer 5 µl of each into a fresh microcentrifuge tube for later flow cytometer analysis. Add 200 µl of 100% (vol/vol) ethanol to the remaining S and G2 phase samples to give samples for DNA extraction.
 - (xvii) Add 500 µl of 1× SYTOX Green in 50 mM sodium citrate to the control 'ALL' sample (from Step 2A(xv)) and the 5-µl S and G2 phase sample aliquots (from Step 2A(xvi)). Incubate these samples in the dark for ≥1 h at RT or at 4 °C overnight to give 1× SYTOX Green re-stained cell samples for flow cytometry.
- **PAUSE POINT** Samples for DNA extraction and flow cytometry from Step 2A(xvi and xvii) can be stored at 4 °C for ≤1 month.
- (xviii) Check the purity of FACS-enriched 'S' and 'G2' cell fractions and compare to the 'ALL' fraction by running the 1× SYTOX Green re-stained cell samples on a BD LSRFortessa X-20 flow cytometer (see Fig. 2). Gate events as in Step 2A(x) and acquire a minimum of 5,000 events gated to assess the sorted samples' purity (Fig. 2c).
 - (xix) Pellet cells in the remaining S and G2 samples by centrifugation at 20,000g for 10 min at RT. Remove the supernatant carefully without disturbing the cell pellet. Proceed to Step 3.

(B) Sync-seq ● **Timing 1.5 d (hands-on time: 3 h)**

- (i) Measure cell density by spectrophotometry or a hemocytometer and inoculate 10^7 cells from the overnight culture into 25 ml of YPAD medium in a 100-ml conical flask; allow to grow to $2\text{--}6 \times 10^7$ cells/ml at 23 °C in a shaking incubator with 200 rpm shaking (~8 h). Check cell density by spectrophotometry or a hemocytometer.
 - (ii) Inoculate 1.5×10^5 cells/ml into 200 ml of YPAD medium in a 1-l conical flask. Allow to grow to 1×10^7 cells/ml at 23 °C in a shaking incubator with 200 rpm shaking (~17 h).
 - (iii) Prepare thirteen 2-ml aliquots of 0.2 M EDTA, 0.1% (wt/vol) sodium azide in 50-ml centrifuge tubes and place tubes tilted at –20 °C overnight.
 - (iv) Divide 10-µl aliquots of 10% (wt/vol) sodium azide into 28 microcentrifuge tubes and place on ice.
 - (v) At 1×10^7 cells/ml, transfer 1 ml of culture into a microcentrifuge tube containing 10 µl of 10% (wt/vol) sodium azide, briefly vortex and place on ice. This is the asynchronous flow cytometry sample.
 - (vi) Add 33 µl of 3 mM alpha factor (final concentration: 0.5 µM) to the culture from Step 2B(ii) and start a timer. (Use a final concentration of 0.2 µM for *BARI* mutant strains.) Repeat every 30 min for 2.5 h (not necessary for *BARI* mutant strains).
- ▲ **CRITICAL STEP** Strains with altered growth rates will need optimization of alpha factor addition times.
- (vii) After 2.5 h, remove 10 µl of culture and inspect cells by light microscopy to monitor cell cycle arrest by counting the fraction of budded cells.
 - (viii) When >95% of cells are unbudded, take the 0-min flow cytometry time point by transferring 1 ml of culture into a microcentrifuge tube containing 10 µl of 10% (wt/vol) sodium azide, vortex and place on ice.

● **TROUBLESHOOTING**

- (ix) To take the 0-min (non-replicating control sample) HTS sample, relocate a pre-prepared 2-ml EDTA/sodium azide aliquot (Step 2B(iii)) from –20 °C to ice and add 100 µl of 10% (wt/vol) sodium azide and 10 ml of the yeast culture to the tube. Immediately shake the

tube vigorously until the frozen pellet dislodges, and then keep on ice until the final HTS sample has been collected.

- (x) Add 2 ml of 20 mg/ml pronase to the culture and start a timer immediately.
- (xi) Take flow cytometry samples (as Step 2B(viii)) every 2.5 min from 5 to 55 min and then at 60, 70, 80, 90, 120 and 180 min after pronase addition.
- (xii) Take HTS samples (as Step 2B(ix)) every 5 min from 5 to 55 min and 90 min after pronase addition.
- (xiii) After taking the 90-min time point, pellet cells from HTS samples by centrifugation at 2,000g for 5 min at RT.
- (xiv) Resuspend cells in 25 ml of Milli-Q water and pellet by centrifugation at 2,000g for 5 min at RT.
- (xv) Remove the supernatant, resuspend the pellet in 1 ml of Milli-Q water and relocate to a screw-capped microcentrifuge tube.
- (xvi) Pellet cells by centrifugation at 13,000g for 30 s at RT, remove the supernatant and store tubes at -20°C .
- (xvii) After the final (180-min) flow cytometry time point, pellet cells from flow cytometry samples by centrifugation at 13,000g for 30 s at RT, remove the supernatant and wash cells in 1 ml of Milli-Q water.
- (xviii) Pellet cells as before (13,000g for 30 s at RT), remove the supernatant and resuspend cells in 1 ml of 70% (vol/vol) ethanol.

■ PAUSE POINT HTS samples can be stored at -20°C for ≤ 1 year. Flow cytometry samples can be stored at 4°C indefinitely.

- (xix) Pellet cells from flow cytometry samples from Step 2B(v, viii and xi) by centrifugation at 13,000g for 1 min at RT and discard the supernatant. Resuspend cells in 1 ml of freshly filter-sterilized 50 mM sodium citrate and pellet by centrifugation at 13,000g for 1 min at RT. Repeat this wash, pellet by centrifugation as before and resuspend cells in 1 ml of 50 mM sodium citrate.
- (xx) Add 25 μl of 10 mg/ml RNaseA. Incubate for 1 h at 37°C in a water bath.
- (xxi) Add 20 μl of 20 mg/ml proteinase K and incubate for 1 h at 55°C in a water bath.
- (xxii) Pellet cells by centrifugation at 13,000g for 1 min at RT. Resuspend cells in 500 μl of $1\times$ SYTOX Green, 0.1% (wt/vol) sodium azide in 50 mM sodium citrate. Incubate for ≥ 1 h at RT or store at 4°C overnight, in the dark.

■ PAUSE POINT SYTOX Green-stained cells can be stored at 4°C for ≤ 1 year.

- (xxiii) Run flow cytometry samples on a BD LSRFortessa X-20 flow cytometer (or another comparable flow cytometer). Gate out doublets based on the forward scatter information as in Fig. 2a. Acquire signal excitation intensity from the 488-nm laser to measure DNA content after SYTOX Green staining.

▲ CRITICAL STEP High-quality data are key to calculating bulk replication in the following step. Acquire 100,000 gated events to assess cell cycle progression of each sample.

- (xxiv) For each time point, extract the median DNA content intensity and normalize the values to have a minimum of 0% (from the 0-min time point) and maximum of 100% (the late time points) DNA content. Plot the normalized values against the time after release that the samples were taken; this shows the increase in bulk genome content and should fit a sigmoid curve (Fig. 3). (Our online tool may be used for the steps described above: <https://dnareplab.shinyapps.io/fitsigmoid/>.) Select time points for HTS based on their genome content (e.g., 10–25% bulk genome content for early S phase, 40–60% for mid S phase and 75–90% for late S phase) (Fig. 3). The bulk replication percentage of the time point samples taken forward for sequencing will be required in Step 61A(viii) or Step 61B(ii).
- (xxv) Proceed to Step 3 with the selected HTS samples (Step 2B(xvi)), as chosen based on bulk genome content and the non-replicating control.

(C) MFA-seq ● Timing 3 d (hands-on time: 2 h)

- (i) Prepare two 2-ml aliquots of 0.2 M EDTA and 0.1% (wt/vol) sodium azide in 50-ml centrifuge tubes and place the tubes tilted at -20°C overnight.
- (ii) Measure the cell density of a 5-ml overnight culture using spectrophotometry or a hemocytometer. Inoculate 6×10^6 cells from the overnight culture into 25 ml of YPAD medium in a 100-ml conical flask and allow to grow to 1.0×10^7 cells/ml at 30°C in a shaking incubator with 200 rpm shaking (~ 5 h).




Fig. 3 | Sync-seq timecourse. Flow cytometry analysis of sync-seq time points. Samples were prepared as described in Step 2B(xix–xxiv). **a**, A typical increase in DNA content in a *Saccharomyces cerevisiae* population released synchronously from alpha factor arrest. **b**, The median DNA content of each time point was plotted against the time the sample was taken, and a sigmoidal curve was fit. Time points for HTS are chosen based on the bulk genome replication: early S phase was captured between 35 and 40 min, mid-S phase between 45 and 50 min and late S phase between 55 and 60 min. This quantification of bulk genome replication provides the normalization values for Steps 61A(viii) and 61B(ii).

- (iii) At 1.0×10^7 cells/ml, remove 1 ml of culture for flow cytometry (process this sample as described in Step 2B(xvii–xxiii)). Relocate a 2-ml EDTA/sodium azide aliquot from -20°C to ice, add 100 μl of 10% (wt/vol) sodium azide and transfer 10 ml of culture into the tube. Immediately shake the tube vigorously until the frozen pellet dislodges, and then place the tube onto ice. This is the replicating sample.
- (iv) Pellet cells from the replicating sample by centrifugation at 2,000g for 5 min at RT.
- (v) Remove the supernatant, resuspend cells in 25 ml of Milli-Q water and pellet by centrifugation at 2,000g for 5 min at RT.
- (vi) Remove the supernatant, resuspend the pellet in 1 ml of Milli-Q water and relocate to a screw-capped microcentrifuge tube.
- (vii) Pellet cells by centrifugation at 13,000g for 30 s at RT, remove the supernatant and store tube at -20°C until the non-replicating sample has been collected and both are ready to proceed to Step 3.
- (viii) Allow the remaining overnight culture to reach stationary phase at 30°C with 200 rpm shaking (~ 48 h).
- (ix) Measure cell density of the overnight culture by spectrophotometry or a hemocytometer. When the cell number stops increasing, relocate a 2-ml EDTA/sodium azide aliquot from -20°C to ice, add 50 μl of 10% (wt/vol) sodium azide and transfer the cell culture into the tube. Immediately shake the tube vigorously until the frozen pellet dislodges, and then place the tube onto ice. This is the non-replicating sample. Process as for the replicating sample, Step 2C(iv–vii).

■ **PAUSE POINT** Samples can be stored at -20°C for ≤ 1 year.

Genomic DNA extraction ● **Timing 3.5 h (hands-on time: 1.5 h)**

- 3 Thaw cell pellets (if required), resuspend cells in 500 μl of solution 1 and add 50 μl of 10 mg/ml zymolyase T20 and 0.5 μl β -mercaptoethanol.
- 4 Incubate samples at 37°C for 30 min in a water bath.

- 5 For sort-seq samples (from Step 2A(xix)), follow option A; for sync-seq samples (from Step 2B (xvi)) or MFA-seq samples (from Step 2C(vii and ix)), follow option B.

(A) **Sort-seq samples**

- (i) Set up the following reactions for each of the samples, incubate at 55 °C for 2 h in a water bath and then proceed to Step 6.

Component	Amount (μl)	Final concentration
Sort-seq sample from Step 4	555.5	–
NaCl (3 M)	17	85 mM
SDS (20% wt/vol)	12.5	0.42% (wt/vol)
Proteinase K (20 mg/ml)	10	0.33 mg/ml
RNase A (10 mg/ml)	5	83.3 μg/ml
Total	600	–

(B) **Sync-seq and MFA-seq samples**

- (i) Pellet cells by centrifugation at 16,000g for 2 min at RT, remove the supernatant and set up the following reactions for each sample. Incubate at 55 °C in a water bath for 2–6 h and then proceed to Step 6.

Component	Amount (μl)	Final concentration
Sync-seq or MFA sample	25	–
Sorbitol (1 M)	50	83.3 mM
Solution 2	500	–
Proteinase K (20 mg/ml)	10	0.33 mg/ml
RNase A (10 mg/ml)	5	83.3 μg/ml
Total	600	–

- 6 Add 0.5 ml of phenol:chloroform:isoamyl alcohol (25:24:1 (vol/vol)) and invert vigorously to mix.
- 7 Centrifuge samples at 16,000g for 5 min at RT and transfer the upper aqueous layer into a fresh microcentrifuge tube.
▲ **CRITICAL** Phenol is toxic and corrosive. Wear PPE and work in a fume hood.
- 8 Add another 0.5 ml of phenol:chloroform:isoamyl alcohol (25:24:1 (vol/vol)), mix well by inversion, centrifuge samples at 16,000g for 5 min at RT and transfer the upper phase into a fresh microcentrifuge tube.
- 9 Add 0.5 ml of chloroform, mix well by inversion, centrifuge samples at 16,000g for 5 min at RT and transfer the upper phase into a fresh microcentrifuge tube.
▲ **CRITICAL STEP** Ensure complete removal of residual phenol. If required, perform a second chloroform extraction step.
- 10 Add 1 ml of 100% (vol/vol) ethanol and invert the tube 10 times. Incubate sort-seq samples at –20 °C for 1 h (this incubation is unnecessary for sync-seq and MFA-seq samples).
- 11 Centrifuge samples at 16,000g for 10 min at RT and remove the supernatant carefully without disturbing the DNA pellet.
- 12 Add 200 μl of 70% (vol/vol) ethanol and invert the tube carefully.
- 13 Remove the supernatant carefully without disturbing the DNA pellet.
- 14 Pulse-spin the sample and remove the remaining ethanol.
- 15 Leave the tube open at RT for 2 min to allow the DNA pellet to dry.
- 16 Resuspend the DNA pellet in 100 μl of TE buffer.
- 17 Measure DNA concentration using Qubit dsDNA HS assay. The yield should be 50–100 ng of DNA for sort-seq samples and 1–10 μg of DNA for sync-seq or MFA-seq samples. Assess DNA purity by NanoDrop spectrophotometer. Pure DNA samples should have absorbance ratios of ~1.8 for 260/280 nm and 2.0–2.2 for 260/230 nm.

? **TROUBLESHOOTING**

- **PAUSE POINT** Store DNA at 4 °C for ≤1 month or at –20 °C for longer storage.




Fig. 4 | Library preparation. Example data from TapeStation fragment analysis for key steps during the library-making procedure (Steps 24–46). Lane names in **a** indicate the stage after which the sample was taken for TapeStation analysis. Intensities are scaled to normalize maxima across the samples; the purple and green lines indicate the upper (1.5 kb) and lower (25 bp) markers, respectively. Note primer bands in adaptor ligation and PCR lanes (marked “*”), which are removed by cleanup steps. The intensity plot in **b** shows the size distribution for a representative final library (as second cleanup lane in **a**), and highlighted in green is the region for which average size is calculated and used for quantification purposes (Step 51). FU, fluorescence units.

Fragmentation of DNA by sonication ● Timing 1.5 h (hands-on time: 45 min)

▲ **CRITICAL** For sequencing using the NextSeq 500 (Illumina) platform, DNA should be fragmented to 300–550 bp. If using other sequencing platforms, the recommended insert size should be checked before fragmentation. We fragment DNA using sonication (Bioruptor Standard); however, fragmentation of DNA could also be achieved using other types of mechanical or enzymatic fragmentation.

▲ **CRITICAL** If other fragmentation methods are used, the conditions for the fragmentation procedure must be optimized to generate 300–550-bp fragments.

18 Decide on the input amount of DNA (5 ng to 1 µg) for library-making. It is best to sonicate more than the minimum amount as some DNA may be lost during sonication and precipitation.

19 Dilute your chosen amount of DNA into TE buffer in a total volume of 400 µl, in a 1.5-ml microcentrifuge tube. We have not encountered problems when using DNA concentrations lower than recommended by the manufacturer of the sonicator.

▲ **CRITICAL STEP** The volume is critical to ensure homogeneous shearing. Take care to dilute samples in TE buffer and not water.

20 Insert sample tubes into the microcentrifuge tube holder of a Bioruptor standard, filling empty holes with blank tubes containing 400 µl of liquid.

21 Sonicate for 20 min, on setting ‘L’, 30 s on, 30 s off, in ice cold water.

22 Precipitate DNA by adding 2× volumes of 100% (vol/vol) ethanol and mix by inversion. Centrifuge for 10–30 min at 13,000g and 4 °C.

23 Carefully discard the supernatant, wash the pellet in 500 µl of 80% (vol/vol) ethanol and air dry. Resuspend the pellet in 29 µl of 1× TE buffer.

24 Check that fragmentation is predominantly between 300 and 550 bp using 2 µl of your sample with a TapeStation (see Fig. 4a; alternatively use Bioanalyzer, Agilent or a PAGE gel), and determine the concentration by Qubit (or another high-sensitivity DNA concentration measurement).

? TROUBLESHOOTING

■ **PAUSE POINT** DNA can be stored overnight at 4 °C. We do not recommend leaving fragmented DNA for longer than this.

End Preparation ● Timing 1.5 h (hands-on time: 10 min)

25 If you have multiple samples, the amount of DNA carried forward into library preparation is determined by the sample with the lowest yield. Take maximal, but equal, amounts of DNA for each sample by diluting samples as appropriate in TE buffer.

- 26 For each sample, add the following to a PCR tube or strip and then mix and pulse spin:

Component	Amount (μl)	Final concentration
500 pg to 1 μg fragmented DNA in TE buffer	25	–
NEBNext Ultra II End Prep Enzyme Mix	1.5	0.05×
NEBNext Ultra II End Prep Reaction Buffer	3.5	0.12×
Total	30	–

- 27 In a thermal cycler with heated lid at 100 °C, incubate the PCR tubes for the following amount of time:

20 °C	30 min
65 °C	30 min
4 °C	hold

▲ **CRITICAL STEP** Continue to the next step ideally straight away, or within 30 min.

Adaptor ligation ● **Timing 1 h (hands-on time: 20 min)**

- 28 Depending on the amount of DNA carried forward into library preparation (Step 25), prepare a dilution (1.25 μl per sample) of NEBNext Adaptor for Illumina in 10 mM Tris-HCl pH 8.0, 10 mM NaCl. Use the following table as a guide:

Amount of DNA	Dilution of adaptor
101 ng to 1 μg	No dilution
5–100 ng	1:10
<5 ng	1:25

- 29 Add the following to each End Prep tube from Step 27 and then mix and pulse spin:

Component	Amount (μl)	Final concentration
End prep sample (Step 27)	30	–
NEBNext Ultra II Ligation Master Mix	15	0.32×
NEBNext Ligation Enhancer	0.5	0.01×
NEBNext Adaptor for Illumina dilution from Step 28	1.25	0.027×
Total	46.75	–

- 30 In a thermal cycler with the heated lid off, incubate tubes for 25 min at 20 °C.
 31 Add 1.5 μl of USER enzyme (NEBNext Multiplex Kit) to each PCR tube and mix.
 32 Incubate 15 min at 37 °C in a thermal cycler with the heated lid at 50 °C.

■ **PAUSE POINT** Samples can be stored overnight at –20 °C; however, we prefer to pause after bead cleanup (Step 40).

Cleanup ● **Timing 1.5–2.5 h (hands-on time: 40 min)**

- 33 Equilibrate AMPure XP beads to RT for 30 min. Vortex for 30 s.
 34 Prepare 400 μl per sample of fresh 80% (vol/vol) ethanol.
 35 Add 43.5 μl (0.9×) of beads to USER reactions from Step 32 and transfer to a 1.5-ml microcentrifuge tube. Vortex the tubes and incubate on a rotating wheel at RT for 10–60 min.
 36 Pulse-spin the tubes and place on a magnetic rack for 5 min until the beads have separated. Discard the supernatant without disturbing the beads.
 37 Pipette 200 μl of freshly made 80% (vol/vol) ethanol into the tube on the side away from the beads. Incubate for 30 s and remove the ethanol. Repeat this 80% (vol/vol) ethanol wash.

- 38 Close the tubes and pulse-spin. Pipette off any remaining ethanol. Air dry for ~30 s or until the beads are still glossy but starting to dry at the edges.
- ▲ **CRITICAL STEP** Both excess ethanol and overly dried beads can lead to poor DNA recovery. If getting poor DNA recovery, we recommend that the air dry time is optimized.
- 39 Elute DNA by adding 8 µl of 10 mM Tris HCl pH 8.0 and vortex well; place in a shaking heat block at 37 °C for 10 min.
- 40 Pulse-spin tubes and place on a magnetic rack for 5 min or until separated (the beads may not separate quite as well as in Step 36). Transfer 7.5 µl to a new PCR tube.
- **PAUSE POINT** Samples can be left overnight at 4 °C.

PCR enrichment of adaptor-ligated DNA ● Timing 1 h (hands-on time: 20 min)

- 41 Decide on a suitable number of PCR cycles for your starting amount of DNA. We use three cycles (the minimum number of cycles) for 450 ng, four cycles for 230 ng and seven cycles for 45 ng.
- 42 Decide on the index primers to use for your samples, taking account of NEB-recommended combinations of primers for pools of six samples or fewer.
- 43 To the PCR tubes with your samples from Step 40, add the following and then mix and pulse spin:

Component	Amount (µl)	Final concentration
Adaptor-ligated samples from Step 40	7.5	-
NEBNext Ultra II Q5 Master mix	12.5	0.5×
Index Primer (NEBNext Multiplex Kit)	2.5	0.1×
Universal primer (NEBNext Multiplex Kit)	2.5	0.1×
Total	25	-

- 44 In a thermal cycler with heated lid at 100 °C, amplify with the following PCR conditions:

Cycle number	PCR conditions	
1	Denature	98 °C, 30 s
2–4*	Denature	98 °C, 10 s;
	anneal/extension	65 °C, 75 s
5	Extension	65 °C, 5 min
6	Hold	4 °C

* Repeat 3+ cycles depending on the amount of input DNA (see Step 41).

Cleanup ● Timing 1.5–2.5 h (hands-on time: 40 min)

- 45 Repeat Steps 33–40 except add 22.5 µl (0.9×) of beads to the reactions from Step 44 and use 33 µl of 0.1× TE buffer for elution. Transfer the final library to a new tube.

Quality-control checks ● Timing 3 h (hands-on time: 1 h)

- 46 Check the read fragment profile of the final library by TapeStation using 2 µl of the final library (see Fig. 4). Alternatively, a Bioanalyzer or PAGE gel can be used.

? TROUBLESHOOTING

- 47 Check the concentration of correctly adaptor-ligated fragments for each library with NEBNext Library Quant Kit for Illumina using a quantitative PCR instrument and appropriate tubes. We use a Qiagene Rotor Gene and PCR tubes, using 20-µl reactions, without ROX, as recommended by the NEBNext Library Quantification protocol. We routinely check 1:100,000 dilutions of final libraries in triplicate; however, initially we recommend running both 1:10,000 and 1:100,000 dilutions of final libraries in triplicate. These dilutions should be made in 1× NEBNext Library Quant Dilution Buffer using a series of 1:100 and 1:10 dilutions.
- 48 Prepare sufficient NEBNext Library Quant Master Mix (with primers) for 15 control reactions, 3 reactions per library and one spare reaction, where each reaction requires 16 µl of the mix. The mix is made by adding a 1/16 volume of NEBNext Library Quant Primer Mix to a 15/16 volume of NEBNext Library Quant Master Mix. For example, for eight libraries prepare 640 µl of NEBNext

Library Quant Master Mix (with primers) by adding 40 µl of NEBNext Library Quant Primer Mix to 600 µl of NEBNext Library Quant Master Mix.

- 49 For a no-template control (1× Library Dilution Buffer), each NEBNext Library Quant Kit DNA standard (10, 1, 0.1 and 0.01 pM) and each library sample, prepare qPCR reactions in triplicate:

Component	Amount (µl)
NEBNext Library Quant Master Mix (with primers) from Step 48	16
No-template control*, DNA standard or library dilution	4
Total	20

* The threshold cycle (Ct) value from the no-template control allows detection of reagent contamination, and it is not used in subsequent quantitation analysis.

- 50 Run the qPCR program and acquire signal:

Cycle number	qPCR conditions
1	Denature 95 °C, 60 s
2–36	Denature 95 °C, 15 s; anneal/extension 63 °C, 45 s acquisition (470 nm/510 nm)

- 51 For each DNA standard, determine the average Ct value from the triplicates. Use these values to plot a standard curve of average Ct against log(pM) and fit a linear regression from which the gradient (m) and intercept with the y-axis (c) can be determined. Calculate the pM concentration for all library dilution triplicates from the Ct value and standard curve gradient (m) and intercept (c) values, using the formula:

$$10^{((c-Ct\text{ value})/m)}$$

Average across triplicates and adjust for the average fragment size for the library (from TapeStation) using the formula:

$$(Triplicate\ average \times 399) / average\ library\ fragment\ size\ (from\ Step\ 46)$$

- 52 Multiply the corrected dilution concentration by the dilution factor (e.g., × 10,000 for a 1:10,000 dilution) and average the concentration across both the 1:10,000 and 1:100,000 dilutions if both dilutions were performed. This gives the final concentration of the correctly adaptor-ligated fragments of the library.


? TROUBLESHOOTING

Steps 51 and 52 can be carried out using NEB's NGS Library Quant web tool: <http://nebiocalculator.neb.com/#!/qPCRlibQnt>.

■ **PAUSE POINT** Libraries can be stored at 4 °C for 2 weeks or at −20 °C for 12 months.

Preparation for, and loading of, the sequencing run ● **Timing 13 h (hands on time: 1 h)**

- 53 Set up a NextSeq run using Basespace (Illumina: <https://basespace.illumina.com/home/index>). Fill in sample names, barcode numbers and cycle number settings. By default, we use single-end sequencing, but the downstream computational pipeline can also be used with paired-end sequencing data.
- 54 The evening before the NextSeq run: follow instructions provided with NextSeq 500/550 High Output Kit v2.5 for cartridge preparation.
- 55 Pool libraries to obtain >5 µl at 4 nM. Optionally, verify the pooled library concentration using the NEBNext Library Quant Kit as described in Steps 47–52.
- 56 Denature, dilute and load the library as recommended for NextSeq500/550 by Illumina. A 75-cycle run takes ~13 h from loading the library to downloading data from Basespace.
- ▲ **CRITICAL STEP** The optimal concentration of library to load for sequencing using a NextSeq 500/550 should be optimized by individual laboratories. We have found that loading a 2.2 pM



pooled library generates a high read number without compromising read quality. We expect to obtain 400–500 million reads where >85% of reads have a quality score $\geq Q30$ per NextSeq run.

■ PAUSE POINT Undiluted pooled libraries can be stored at 4 °C for 2 weeks or at –20 °C for 12 months.

◀ **Fig. 5 | Anticipated results.** Anticipated results and graphs generated at key steps through the Repliscope R package. **a**, Typical reads per bin distribution of a non-replicating sample (Step 61A(iii) or 61B(ii)), after removing values that are smaller than $0.25 \times$ median. The median is shown as a pink line across the panels. Most bins have similar values, and, therefore, bins containing fewer reads are detected as outliers (highlighted in gray) using the interquartile range method. Usually, there is no need to remove them. **b**, Typical reads per bin distribution of a replicating sample (Step 61A(v) or 61B(ii)) after applying $0.25 \times$ median filtering. A wider distribution of reads than in non-replicating samples is expected. **c**, Chromosome (chr) VIII sequencing coverage of the sample in **a** (Step 61A(iv) or 61B(ii)). Data clustering close to the genome-wide median (pink line) is indicative of a good non-replicating sample. **d**, chrVIII sequencing coverage of the sample in **b** (Step 61A(v) or 61B(ii)). Peaks and valleys are indicative of a good replicative sample. **e**, Sequencing depth-adjusted ratio distribution centers around 1 (Step 61A(vii) or 61B(ii)). Region 1–2 is highlighted in green. **f**, A typical sort-seq sample with ratio values auto-normalized to fit a biologically meaningful scale (Step 61A(viii) or 61B(ii)). Region 1–2 is highlighted in green. **g**, Sort-seq replication profile of *Saccharomyces cerevisiae* chrVIII (blue points) with fitted spline smoothing (gray line). **h**, Sync-seq replication profile of *Saccharomyces cerevisiae* chrXI. **i**, MFA-seq replication profile of *H. volcanii* chromosome. g–i: Step 61A (ix and x) or 61B(ii).

Sequencing read alignment and generation of raw read counts ● Timing 2–8 h (hands-on time: 30 min)

57 Login to your Illumina account. Check run quality information and download the fastq sequencing files. Alternatively, follow instructions provided by your sequencing facility to download fastq sequencing files.

? TROUBLESHOOTING

58 Use FastQC to make sure that the sequencing files are intact and to check the quality of the data generated from the sequencing run. See Table 1 for more information.

? TROUBLESHOOTING

59 If using a reference genome for the first time, it needs to be indexed with Bowtie 2. Obtain a multiple sequence fasta file containing the genome contigs (e.g., *sacCer3* is available from https://www.ncbi.nlm.nih.gov/assembly/GCF_000146045.2/) and use the following command to build a Bowtie 2 index:

```
$ bowtie2-build -f /path/to/genome.fasta genomeName
```

This will create indexed genome files with bw2 or bw2l extensions in the directory where the command was executed.

60 To align the reads to the index genome (created in Step 59) and count the reads in genomic windows (bins), download the localMapper script from GitHub (<https://github.com/DNAReplicationLab/localMapper>) and make executable:

```
$ sudo chmod +x filePath/localMapper.sh
```

Navigate to the directory containing the sequencing files and follow option A for single-end sequencing and Option B for paired-end sequencing. Use -h argument to print more information about available options. By default, the script will create a subdirectory with the provided sample name. Once the script has finished, the sample_name/processed directory will contain a non-empty bed file with read counts in genomic bins of the size specified by the -w argument. We typically use 1-kb bins for small genomes (e.g., yeasts) and 50-kb bins for larger genomes (e.g., mammals); this is discussed further in the Experimental design.

(A) Single-end sequencing

(i) Run the script as follows:

```
$ filePath/localMapper.sh -g <genome name> -U <fastq file(s)> -s
<sample name> -c <number of threads> -m <memory per thread> -w
<window size in bp>
```

(B) Paired-end sequencing

(i) Run the script as follows:

```
$ filePath/localMapper.sh -g <genome name> -1 <first mate fastq file
(s)> -2 <second mate fastq file(s)> -s <sample name> -c <number of
threads> -m <memory per thread> -w <window size in bp>
```

Generate, normalize and plot the ratios ● **Timing 5–30 min**

61 To obtain replication profiles, the count files from replicating samples have to be normalized using counts from non-replicating samples to produce copy-number ratios. The ratios are normalized based on the sequencing depth, followed by a dynamic range adjustment. The adjusted ratios are visualized and may be compared using statistics. All these steps can be performed using the Repliscope package available from CRAN and GitHub (<https://github.com/DNAReplicationLab/Repliscope>). To perform the analysis via the command line, follow option A. To perform the analysis interactively in a web browser, follow option B.

(A) Data analysis using command line functions

(i) Launch R and load the Repliscope package:

```
$ R> library(Repliscope)
```

(ii) Load the bed file from the non-replicative sample produced in Step 60.

```
> nonRepSample <- loadBed('filePath/sample.bed')
```

(iii) Visualize the loaded read counts as a boxplot to check for aneuploidy and outliers. All chromosomes, except for mitochondria, should cluster around the genome-wide median, as shown in Fig. 5a. Aneuploidy can be identified if one or more chromosomes are higher or lower than the median line (Fig. 6a). A group of outliers on a single chromosome above the median indicates a segmental duplication (Fig. 6b).

```
> plotBed(nonRepSample)
```

If outliers are detected, use the `rmOutliers()` function to remove them. As a standard, we remove bins that contain less than a quarter of the genomic median counts value; these tend to produce noisier ratios:

```
> nonRepSample <- rmOutliers(nonRepSample, 'median', loLim=0.25)
```

? TROUBLESHOOTING

(iv) Visualize read counts as a scatterplot using the following command:

```
> plotCoverage(nonRepSample)
```

This provides a visual check that read count values are even across chromosomes (Fig. 5c). Curvature on the plot may indicate that there was ongoing replication in the sample, and it may be necessary to use a different non-replicative control.

? TROUBLESHOOTING

The plot also allows aneuploidy to be identified if one or more chromosomes are higher or lower than the median line (Fig. 6c), while a segmental duplication is seen as a region of a chromosome that is abruptly higher than the median (Fig. 6d).

(v) Repeat Step 61A(ii–iv) using the replicative sample counts file and replace ‘nonRepSample’ with ‘repSample’ in the commands (Fig. 5b,d).

(vi) Create the sequencing depth–normalized count ratio:

```
> theRatio <- makeRatio(repSample, nonRepSample)
```

This creates a dataframe that contains a ‘ratio’ column with calculated values adjusted by total read number.

(vii) Visualize the calculated ratio values as a histogram (Fig. 5e).

```
> plotRatio(theRatio$ratio)
```

The values should be mostly between 0.5 and 1.5 and center around 1. Outliers may be removed using the `trimRatio()` function.




Fig. 6 | Chromosomal aberrations. Two examples of samples showing chromosomal aberrations. Affected chromosomes are labelled in red. **a** and **c**, Evidence of chrXI and chrXII aneuploidy from a sort-seq experiment on a diploid *Saccharomyces cerevisiae* strain. Data are plotted as reads per bin, as a boxplot view (**a**; Step 61A(iii) or 61B(ii)) and a whole chromosome view (**c**; Step 61A(iv) or 61B(ii)). The region with missing data on chrXII is the location of an rDNA repeat array. **b** and **d**, Evidence of segmental duplication of a region of chrIV. Data are plotted as reads per bin, as a boxplot view (**b**) and a whole chromosome view (**d**; Step 61A(v) or 61B(ii)). The pink line across the panels shows the whole genome median count. The gray rectangle in **b** highlights outliers using the interquartile range method.

(viii) Normalize the ratio to fit a biologically relevant scale, either using the default autonormalization feature or by providing a normalization factor manually (Box 3). In case of sync-seq, the normalization factor is calculated from bulk replication values obtained in Step 2B(xxiv) using the formula $1 + [\text{bulk genome replication as a decimal}]$. For example, use the following command to normalize a sync-seq sample with 35% of bulk genome replication:

```
> theRatio <- normaliseRatio(theRatio, rFactor=1.35)
```

In case of sort-seq, the normalization factor is determined algorithmically, by minimizing the sum of values outside of the 1 to 2 scale. To do so, run the function without additional arguments:

```
> theRatio <- normaliseRatio(theRatio)
```

For MFA-seq, the maximum relative copy number can be calculated as $1 + [\text{proportion of cells in S phase}]$. To run the minimizing function with the new maximum (e.g., 1.25), run the function like this:

```
> theRatio <- normaliseRatio(theRatio, upperLimit=1.25)
```

Visualize the normalized ratio using the command from Step 61A(vii) (Fig. 5f).

(ix) Plot the resulting replication profile:

```
> plotGenome(theRatio)
```

? TROUBLESHOOTING

(x) To plot smoothed data, run the spline smoothing function:

```
> theRatio <- smoothRatio(theRatio)
```

Repeat Step 61A(ix) to plot smoothed ratio values (Fig. 5g).

Box 4 | Data interpretation and custom analysis

Following our protocol, users will generate replication timing profiles for their cell type(s) of interest, and can identify genomic regions with a statistically significant difference between two samples. To facilitate comparison of sync-seq datasets, median replication time (or Trep) can be calculated from a series of time points throughout S phase, as described previously^{4,9}. We provide this functionality within the Repliscope package. However, depending on the specific research question, there may be other custom analyses that the user wishes to perform, some of which we suggest here.

Replication profiles allow localization of active replication origins. These are evident as peaks in replication profiles and therefore can be identified with peak calling tools, such as in ref.⁴. Motif searches can be performed as described previously⁴, to look for sequence motifs associated with replication timing profile peaks and, therefore, origins. In addition, the proportion of forks moving in either direction genome-wide can be determined from the gradients of replication timing profiles, and therefore can reveal firing efficiencies for replication origins^{9,59}. Evolutionary comparisons can be made between replication timing profiles from related species of various evolutionary distance (e.g., using the liftOver tool or considering regions of various levels of conservation)^{12,16,21}. Replication timing profiles can also be compared to genomic variation, such as single-nucleotide polymorphisms (SNPs) and alleles identified that associate with replication timing differences¹⁷. Of note, replication timing profiles can generate genotype information such as large-scale rearrangements (see Box 2 and ref.¹⁴). Finally, as a genome-wide dataset, replication timing profiles can be compared with other genome-wide datasets of interest, such as ChIP-seq, RNA-seq or Hi-C.

- (xi) To compare two replication profiles, each generated as described above in Steps 60–61A (viii) and named ‘theRatio1’ and ‘theRatio2’, use the following command:

```
> comparedRatios <- compareRatios(theRatio1,theRatio2)
```

The resulting dataframe will contain both replication profiles and an additional ‘p.value’ column. It can be plotted as in Step 61A(ix).

- (xii) To plot multiple replication profiles (each generated as described above in Steps 60–61A (viii) and named ‘ratio1’, ‘ratio2’, ‘ratio3’, etc.) on the same plot (e.g., for sync-seq), combine them first into a single dataframe:

```
> combinedReplicationProfiles <- rbind(ratio1,ratio2, ratio3...)
```

These replication profiles can be plotted using the command in Step 61A(ix). An example plot is shown in Fig. 5h.

- (B) Data analysis using our Shiny app

- (i) Launch R and load and launch Repliscope in the interactive mode:

```
$ R> library(Repliscope)> runGUI()
```

- (ii) Follow instructions on the webpages of the app to create, normalize and plot the replicating to non-replicating sample counts ratios (Supplementary Note).

? TROUBLESHOOTING

- 62 Depending on the specific application users may then wish to perform custom analyses on the data generated, some of which we suggest in Box 4.

Troubleshooting

Troubleshooting advice can be found in Table 1.

Table 1 | Troubleshooting table

Step	Problem	Possible reason	Solution
2A(x)	Cell flocculation Cell debris	Insufficient cell separation	Perform additional sonication (Step 2A(ix))
		Excessive sonication	Repeat sample preparation and separate cells by vortexing
	Insufficient discrimination of cell cycle stages by FACS	Bacterial infection in yeast culture	Repeat sample preparation, ensure sterility of reagents and consider adding an antibacterial agent
		Incomplete proteinase K or RNase A digest	Check to ensure use of correct proteinase K and RNase A digestion conditions
		Insufficient SYTOX Green labeling	Check to ensure use of correct SYTOX Green labeling conditions
		Misaligned FACS lasers	Check alignment of FACS lasers

Table continued

Table 1 (continued)

Step	Problem	Possible reason	Solution
2B (viii)	Incomplete alpha factor arrest	Non-optimal arrest conditions for strain	Check <i>BAR1</i> status of yeast. Optimize alpha factor arrest for your strain of interest
17	RNA contamination	Incomplete digestion by RNase A	Check to ensure use of correct RNaseA digestion conditions; if necessary, additional RNase A digestions can be performed (Step 5)
	Low amount of extracted DNA	Incomplete cell recovery during centrifugation Lost pellets Incomplete DNA resuspension	Centrifuge at recommended speeds Take care at Step 2A(xix) and Steps 11, 13 and 14 Leave overnight in the refrigerator or heat at 60 °C (Step 16)
24	Low 260/230-nm ratio Lower-than-expected sized fragments	Contaminants in DNA Sonicated in water	Include additional chloroform step Sonicate in TE buffer (Step 19)
24	Larger-than-expected sized fragments	Sonicated for too long or using too high power	Check using correct sonication conditions or optimize if using alternative sonicator (Step 21)
46	Shorter-than-expected library size	Sonicated in smaller volume than 400 µl Adaptor dimer contamination	Sonicate samples in 400 µl; the volume is more critical than the concentration (Step 19) Repeat AMPure XP bead cleanup (Step 45) or repeat library-making steps with higher concentration of DNA and correct adaptor concentration (Step 28)
52	Low library concentration	Starting size of fragmented DNA too small Overestimation of DNA input	See Troubleshooting for Step 24 Check 260/230-nm and 260/280-nm ratios (Step 17) Use higher starting concentration of DNA (Step 18) Increase the number of PCR cycles during library amplification (Steps 41–44)
57	Low read number sequenced but high Q30 score	Low concentration of pooled library loaded onto NextSeq	Optimize concentration of loaded pooled library (Step 55)
57 and 58	High read number but low Q30 score and FastQC quality scores	Too high concentration of pooled library loaded, causing overformation of clusters	Optimize concentration of loaded pooled library (Step 55)
58	High amount of sequence duplication	Low sequencing complexity due to overamplification of limited library	Use higher starting concentration of DNA (Step 18)
	High adaptor content	Adaptor dimer contamination	See Troubleshooting for Step 46
	Presence of index sequences seen in kmer content	Index sequences not removed appropriately	Manually trim reads to remove index sequences
61A(iii) and 61B(ii)	Different average read count for one or more chromosomes	Possible aneuploidy or segmental duplication	See Box 2
61A(iv) and 61B(ii)	Discontinuity in genomic reads	Possible genomic rearrangement	See Box 2
	Increased recovery of chromosome ends in non-replicating sample	Biases from DNA purification, a known problem with some DNA extraction kits	Purify DNA by phenol-chloroform extraction and ethanol precipitation
	Non-replicating sample is not flat across the genome.	The sample is not completely non-replicating and is not suitable to use for normalization.	Use a different non-replicating sample from the same run if possible or repeat sample collection
61A(ix) and 61B(ii)	Low dynamic range (i.e., data points do not reach the expected maxima)	Low dynamic range in original cell sample Intercellular stochasticity of replication	Poor sort purity or gate setting (Step 2A(x–xii)); e.g., contamination of the S phase sample with G2 cells reduces the dynamic range, as does the S phase gate not capturing the full range of S phase cells. Poor cell cycle synchronization: optimize (Step 2B(vi–xi)) Other techniques can be useful to address stochasticity, such as D-NA ⁴³ , which identifies DNA replication on single molecules

Timing

Step 1, growing primary cell culture: 17 h (hands-on time: 10 min)
Step 2, obtaining replicating and non-replicating cells: 1.5–3 d
Step 2A(i–xix), sort-seq: 2 d (hands-on time: 4 h)
Step 2B(i–xxv), sync-seq: 1.5 d (hands-on time: 3 h)
Step C(i–ix), MFA-seq: 3 d (hands-on time: 2 h)
Steps 3–17, genomic DNA extraction: 3.5 h (hands-on time: 1.5 h)
Steps 18–24, fragmentation of DNA by sonication: 1.5 h (hands-on time: 45 min)
Steps 25–27, end preparation: 1.5 h (hands-on time: 10 min)
Steps 28–32, adaptor ligation: 1 h (hands-on time: 20 min)
Steps 33–40, cleanup: 1.5–2.5 h (hands-on time: 40 min)
Steps 41–44, PCR enrichment of adaptor-ligated DNA: 1 h (hands-on time: 20 min)
Step 45, cleanup: 1.5–2.5 h (hands-on time: 40 min)
Steps 46–52, quality-control checks: 3 h (hands-on time: 1 h)
Steps 53–56, preparation for, and loading of, the sequencing run: 13 h (hands on time: 1 h)
Steps 57–60, sequencing read alignment and generation of raw read counts: 2–8 h (hands-on time: 30 min)
Step 61, generate, normalize and plot the ratios: 5–30 min
Step 62, custom analysis: user defined

Anticipated results

Upon successful sequencing and read mapping, non-replicating samples should have a fairly flat profile, featuring a coefficient of variation <20%. This is evidenced by a similar distribution of read number between chromosomes (Fig. 5a) and sequencing coverage within the chromosomes (Fig. 5c). Some windows may contain fewer than expected reads due to filtering out reads mapping to non-unique regions. Replicating samples, on the other hand, should have a wider distribution of read numbers (coefficient of variation around 100%) as evidenced by the increased interquartile range in box plots (Fig. 5b) and curved coverage profiles across chromosomes (Fig. 5d). Once adjusted by the sequencing depth ratio between replicating and non-replicating samples, data points should be distributed around 1 (Fig. 5e). Further normalization will adjust the ratio distribution to fit a biologically relevant scale, usually to fit a 1 to 2 relative copy number scale (Fig. 5f). Normalized ratios, when plotted as a function of genomic position, generate the replication profile for the sample (Fig. 5g–i). The curvature of the profile should be similar to the replicating sample read coverage and inversely proportional to the replication time of the DNA region⁴.

Reporting Summary

Further information on research design is available in the Nature Research Reporting Summary linked to this article.

Data availability

All data generated or analyzed during this study are publicly available from the NCBI GEO functional genomics data repository with the following accession numbers: GSE42243, GSE48212 (Fig. 5); GSE135178 (Fig. 6).

Code availability

The custom bash script required for the analysis, as well as the script to download and analyze example data, are available from GitHub (<https://github.com/DNAReplicationLab/localMapper/>). The R package Repliscope described here is available from CRAN (<https://cran.r-project.org/web/packages/Repliscope/>) and GitHub (<https://github.com/DNAReplicationLab/Repliscope/>). We have also provided an official Ubuntu Desktop 18.04 LTS installation disk image with all the software required for the analysis (<https://ln1.path.ox.ac.uk/groups/nieduszynski/Replibuntu/Replibuntu-18.04.0-amd64.iso.gz>). The code in this manuscript has been peer reviewed.

References

1. Fragkos, M., Ganier, O., Coulombe, P. & Mechali, M. DNA replication origin activation in space and time. *Nat. Rev. Mol. Cell. Biol.* **16**, 360–374 (2015).
2. Newman, T. J., Mamun, M. A., Nieduszynski, C. A. & Blow, J. J. Replisome stall events have shaped the distribution of replication origins in the genomes of yeasts. *Nucleic Acids Res.* **41**, 9705–9718 (2013).
3. Rhind, N. & Gilbert, D. M. DNA replication timing. *Cold Spring Harb. Perspect. Biol.* **5**, a010132 (2013).
4. Müller, C. A. et al. The dynamics of genome replication using deep sequencing. *Nucleic Acids Res.* **42**, e3 (2014).
5. Yabuki, N., Terashima, H. & Kitada, K. Mapping of early firing origins on a replication profile of budding yeast. *Genes Cells* **7**, 781–789 (2002).
6. Raghuraman, M. K. et al. Replication dynamics of the yeast genome. *Science* **294**, 115–121 (2001).
7. Peace, J. M., Villwock, S. K., Zeytounian, J. L., Gan, Y. & Aparicio, O. M. Quantitative BrdU immunoprecipitation method demonstrates that Fkh1 and Fkh2 are rate-limiting activators of replication origins that reprogram replication timing in G1 phase. *Genome Res.* **26**, 365–375 (2016).
8. Sekedat, M. D. et al. GINS motion reveals replication fork progression is remarkably uniform throughout the yeast genome. *Mol. Syst. Biol.* **6**, 353 (2010).
9. Hawkins, M. et al. High-resolution replication profiles define the stochastic nature of genome replication initiation and termination. *Cell Rep.* **5**, 1132–1141 (2013).
10. Claycomb, J. M. & Orr-Weaver, T. L. Developmental gene amplification: insights into DNA replication and gene expression. *Trends Genet.* **21**, 149–162 (2005).
11. Koren, A., Soifer, I. & Barkai, N. MRC1-dependent scaling of the budding yeast DNA replication timing program. *Genome Res.* **20**, 781–790 (2010).
12. Müller, C. A. & Nieduszynski, C. A. Conservation of replication timing reveals global and local regulation of replication origin activity. *Genome Res.* **22**, 1953–1962 (2012).
13. Skovgaard, O., Bak, M., Lobner-Olesen, A. & Tommerup, N. Genome-wide detection of chromosomal rearrangements, indels, and mutations in circular chromosomes by short read sequencing. *Genome Res.* **21**, 1388–1393 (2011).
14. Hawkins, M., Malla, S., Blythe, M. J., Nieduszynski, C. A. & Allers, T. Accelerated growth in the absence of DNA replication origins. *Nature* **503**, 544–547 (2013).
15. Samson, R. Y. et al. Specificity and function of archaeal DNA replication initiator proteins. *Cell Rep.* **3**, 485–496 (2013).
16. Agier, N. et al. The evolution of the temporal program of genome replication. *Nat. Commun.* **9**, 2199 (2018).
17. Koren, A. et al. Genetic variation in human DNA replication timing. *Cell* **159**, 1015–1026 (2014).
18. Batrakou, D. G., Heron, E. D. & Nieduszynski, C. A. Rapid high-resolution measurement of DNA replication timing by droplet digital PCR. *Nucleic Acids Res.* **46**, e112 (2018).
19. Shen, Y. et al. Deep functional analysis of synII, a 770-kilobase synthetic yeast chromosome. *Science* **355**, eaaf4791 <https://doi.org/10.1126/science.aaf4791> (2017).
20. Annaluru, N. et al. Total synthesis of a functional designer eukaryotic chromosome. *Science* **344**, 55–58 (2014).
21. Müller, C. A. & Nieduszynski, C. A. DNA replication timing influences gene expression level. *J. Cell Biol.* **216**, 1907–1914 (2017).
22. Natsume, T. et al. Kinetochore coordinate pericentromeric cohesion and early DNA replication by Cdc7-Dbf4 kinase recruitment. *Mol. Cell* **50**, 661–674 (2013).
23. Kedziora, S. et al. Rif1 acts through Protein Phosphatase 1 but independent of replication timing to suppress telomere extension in budding yeast. *Nucleic Acids Res.* **46**, 3993–4003 (2018).
24. Foss, E. J. et al. SIR2 suppresses replication gaps and genome instability by balancing replication between repetitive and unique sequences. *Proc. Natl Acad. Sci. USA* **114**, 552–557 (2017).
25. Voichkek, Y., Bar-Ziv, R. & Barkai, N. Expression homeostasis during DNA replication. *Science* **351**, 1087–1090 (2016).
26. Rudolph, C. J., Upton, A. L., Stockum, A., Nieduszynski, C. A. & Lloyd, R. G. Avoiding chromosome pathology when replication forks collide. *Nature* **500**, 608–611 (2013).
27. Tanaka, S., Nakato, R., Katou, Y., Shirahige, K. & Araki, H. Origin association of Sld3, Sld7, and Cdc45 proteins is a key step for determination of origin-firing timing. *Curr. Biol.* **21**, 2055–2063 (2011).
28. Green, B. M., Finn, K. J. & Li, J. J. Loss of DNA replication control is a potent inducer of gene amplification. *Science* **329**, 943–946 (2010).
29. Tanny, R. E., MacAlpine, D. M., Blitzblau, H. G. & Bell, S. P. Genome-wide analysis of re-replication reveals inhibitory controls that target multiple stages of replication initiation. *Mol. Biol. Cell* **17**, 2415–2423 (2006).
30. Kramara, J., Osia, B. & Malkova, A. Break-induced replication: the where, the why, and the how. *Trends Genet.* **34**, 518–531 (2018).
31. Alvino, G. M. et al. Replication in hydroxyurea: it's a matter of time. *Mol. Cell. Biol.* **27**, 6396–6406 (2007).
32. Macheret, M. & Halazonetis, T. D. Monitoring early S-phase origin firing and replication fork movement by sequencing nascent DNA from synchronized cells. *Nat. Protoc.* **14**, 51–67 (2019).
33. Macheret, M. & Halazonetis, T. D. Intragenic origins due to short G1 phases underlie oncogene-induced DNA replication stress. *Nature* **555**, 112–116 (2018).

34. Marchal, C. et al. Genome-wide analysis of replication timing by next-generation sequencing with E/L Repli-seq. *Nat. Protoc.* **13**, 819–839 (2018).
35. Siow, C. C., Nieduszynska, S. R., Müller, C. A. & Nieduszynski, C. A. OriDB, the DNA replication origin database updated and extended. *Nucleic Acids Res.* **40**, D682–D686 (2012).
36. Prioleau, M. N. & MacAlpine, D. M. DNA replication origins—where do we begin? *Genes Dev.* **30**, 1683–1697 (2016).
37. Jinks-Robertson, S. & Klein, H. L. Ribonucleotides in DNA: hidden in plain sight. *Nat. Struct. Mol. Biol.* **22**, 176–178 (2015).
38. Petryk, N. et al. Replication landscape of the human genome. *Nat. Commun.* **7**, 10208 (2016).
39. McGuffee, S. R., Smith, D. J. & Whitehouse, I. Quantitative, genome-wide analysis of eukaryotic replication initiation and termination. *Mol. Cell* **50**, 123–135 (2013).
40. Feng, W. et al. Genomic mapping of single-stranded DNA in hydroxyurea-challenged yeasts identifies origins of replication. *Nat. Cell Biol.* **8**, 148–155 (2006).
41. Tiengwe, C. et al. Genome-wide analysis reveals extensive functional interaction between DNA replication initiation and transcription in the genome of *Trypanosoma brucei*. *Cell Rep.* **2**, 185–197 (2012).
42. Donaldson, A. D. & Nieduszynski, C. A. Genome-wide analysis of DNA replication timing in single cells: Yes! We're all individuals. *Genome Biol.* **20**, 111 (2019).
43. Müller, C. A. et al. Capturing the dynamics of genome replication on individual ultra-long nanopore sequence reads. *Nat. Methods* **16**, 429–436 (2019).
44. Haase, S. B. & Reed, S. I. Improved flow cytometric analysis of the budding yeast cell cycle. *Cell Cycle* **1**, 132–136 (2002).
45. Langmead, B. & Salzberg, S. L. Fast gapped-read alignment with Bowtie 2. *Nat. Methods* **9**, 357–359 (2012).
46. Li, H. et al. The Sequence Alignment/Map format and SAMtools. *Bioinformatics* **25**, 2078–2079 (2009).
47. Quinlan, A. R. & Hall, I. M. BEDTools: a flexible suite of utilities for comparing genomic features. *Bioinformatics* **26**, 841–842 (2010).
48. Daigaku, Y. et al. A global profile of replicative polymerase usage. *Nat. Struct. Mol. Biol.* **22**, 192–198 (2015).
49. Dobin, A. et al. STAR: ultrafast universal RNA-seq aligner. *Bioinformatics* **29**, 15–21 (2013).
50. Xu, J. et al. Genome-wide identification and characterization of replication origins by deep sequencing. *Genome Biol.* **13**, R27 (2012).
51. Heichinger, C., Penkett, C. J., Bahler, J. & Nurse, P. Genome-wide characterization of fission yeast DNA replication origins. *EMBO J.* **25**, 5171–5179 (2006).
52. Agier, N., Romano, O. M., Touzain, F., Cosentino Lagomarsino, M. & Fischer, G. The spatiotemporal program of replication in the genome of *Lachancea kluyveri*. *Genome Biol. Evol.* **5**, 370–388 (2013).
53. Di Rienzi, S. C. et al. Maintaining replication origins in the face of genomic change. *Genome Res.* **22**, 1940–1952 (2012).
54. Kim, H. S. Genome-wide function of MCM-BP in *Trypanosoma brucei* DNA replication and transcription. *Nucleic Acids Res.* **47**, 634–647 (2019).
55. Devlin, R. et al. Mapping replication dynamics in *Trypanosoma brucei* reveals a link with telomere transcription and antigenic variation. *Elife* **5**, e12765 <https://doi.org/10.7554/eLife.12765> (2016).
56. Marques, C. A., Dickens, N. J., Paape, D., Campbell, S. J. & McCulloch, R. Genome-wide mapping reveals single-origin chromosome replication in *Leishmania*, a eukaryotic microbe. *Genome Biol.* **16**, 230 (2015).
57. Lundgren, M., Andersson, A., Chen, L., Nilsson, P. & Bernander, R. Three replication origins in *Sulfolobus* species: synchronous initiation of chromosome replication and asynchronous termination. *Proc. Natl Acad. Sci. USA* **101**, 7046–7051 (2004).
58. Woodfine, K. et al. Replication timing of the human genome. *Hum. Mol. Genet.* **13**, 191–202 (2004).
59. de Moura, A. P., Retkute, R., Hawkins, M. & Nieduszynski, C. A. Mathematical modelling of whole chromosome replication. *Nucleic Acids Res.* **38**, 5623–5633 (2010).

Acknowledgements

We thank Amanda Williams and Becky Busby (Zoology Sequencing Facility) for help with the NextSeq 500, and Michal Maj and Line Eriksen (Sir William Dunn School Flow Cytometry Facility) for their help with FACS.

Author contributions

All authors wrote and edited the manuscript.

Competing interests

The authors declare no competing interests.

Additional information

Supplementary information is available for this paper at <https://doi.org/10.1038/s41596-019-0287-7>.

Correspondence and requests for materials should be addressed to C.A.N.

Reprints and permissions information is available at www.nature.com/reprints.

Publisher's note Springer Nature remains neutral with regard to jurisdictional claims in published maps and institutional affiliations.

Received: 27 September 2018; Accepted: 16 December 2019;
Published online: 12 February 2020

Related links

Key references using this protocol

Hawkins, M. et al. *Nature* **503**, 544–547 (2013): <https://doi.org/10.1038/nature12650>

Müller, C. A. et al. *Nucleic Acids Res.* **42**, e3 (2014): <https://doi.org/10.1093/nar/gkt878>

Müller, C. A. & Nieduszynski, C. A. *J. Cell Biol.* **216**, 1907–1914 (2017): <https://doi.org/10.1083/jcb.201701061>

Reporting Summary

Nature Research wishes to improve the reproducibility of the work that we publish. This form provides structure for consistency and transparency in reporting. For further information on Nature Research policies, see [Authors & Referees](#) and the [Editorial Policy Checklist](#).

Statistics

For all statistical analyses, confirm that the following items are present in the figure legend, table legend, main text, or Methods section.

- | | |
|-------------------------------------|--|
| n/a | Confirmed |
| <input checked="" type="checkbox"/> | <input type="checkbox"/> The exact sample size (<i>n</i>) for each experimental group/condition, given as a discrete number and unit of measurement |
| <input checked="" type="checkbox"/> | <input type="checkbox"/> A statement on whether measurements were taken from distinct samples or whether the same sample was measured repeatedly |
| <input type="checkbox"/> | <input checked="" type="checkbox"/> The statistical test(s) used AND whether they are one- or two-sided
<i>Only common tests should be described solely by name; describe more complex techniques in the Methods section.</i> |
| <input checked="" type="checkbox"/> | <input type="checkbox"/> A description of all covariates tested |
| <input type="checkbox"/> | <input checked="" type="checkbox"/> A description of any assumptions or corrections, such as tests of normality and adjustment for multiple comparisons |
| <input type="checkbox"/> | <input checked="" type="checkbox"/> A full description of the statistical parameters including central tendency (e.g. means) or other basic estimates (e.g. regression coefficient) AND variation (e.g. standard deviation) or associated estimates of uncertainty (e.g. confidence intervals) |
| <input checked="" type="checkbox"/> | <input type="checkbox"/> For null hypothesis testing, the test statistic (e.g. <i>F</i> , <i>t</i> , <i>r</i>) with confidence intervals, effect sizes, degrees of freedom and <i>P</i> value noted
<i>Give P values as exact values whenever suitable.</i> |
| <input checked="" type="checkbox"/> | <input type="checkbox"/> For Bayesian analysis, information on the choice of priors and Markov chain Monte Carlo settings |
| <input checked="" type="checkbox"/> | <input type="checkbox"/> For hierarchical and complex designs, identification of the appropriate level for tests and full reporting of outcomes |
| <input checked="" type="checkbox"/> | <input type="checkbox"/> Estimates of effect sizes (e.g. Cohen's <i>d</i> , Pearson's <i>r</i>), indicating how they were calculated |

Our web collection on [statistics for biologists](#) contains articles on many of the points above.

Software and code

Policy information about [availability of computer code](#)

Data collection FlowJo 10.5.3, Basespace.

Data analysis Ubuntu 18.04.01 LTS, FastQC v0.11.5, bowtie2 2.3.4.1, samtools 1.7, bedtools 2.26.0, picard-tools 2.8.1, R 3.4.4, Repliscope 1.0.0, ggplot2 3.0.0, shiny 1.1.0, colourpicker 1.0, Google Chrome 70.0.3538.77

For manuscripts utilizing custom algorithms or software that are central to the research but not yet described in published literature, software must be made available to editors/reviewers. We strongly encourage code deposition in a community repository (e.g. GitHub). See the Nature Research [guidelines for submitting code & software](#) for further information.

Data

Policy information about [availability of data](#)

All manuscripts must include a [data availability statement](#). This statement should provide the following information, where applicable:

- Accession codes, unique identifiers, or web links for publicly available datasets
- A list of figures that have associated raw data
- A description of any restrictions on data availability

All data generated or analysed during this study are included in this published article (and its supplementary information files).

Field-specific reporting

Please select the one below that is the best fit for your research. If you are not sure, read the appropriate sections before making your selection.

- ☒ Life sciences ☐ Behavioural & social sciences ☐ Ecological, evolutionary & environmental sciences

Life sciences study design

All studies must disclose on these points even when the disclosure is negative.

Sample size	N/A
Data exclusions	N/A
Replication	N/A
Randomization	N/A
Blinding	N/A

Reporting for specific materials, systems and methods

We require information from authors about some types of materials, experimental systems and methods used in many studies. Here, indicate whether each material, system or method listed is relevant to your study. If you are not sure if a list item applies to your research, read the appropriate section before selecting a response.

Materials & experimental systems		Methods	
n/a	Involved in the study	n/a	Involved in the study
<input checked="" type="checkbox"/>	<input type="checkbox"/> Antibodies	<input checked="" type="checkbox"/>	<input type="checkbox"/> ChIP-seq
<input checked="" type="checkbox"/>	<input type="checkbox"/> Eukaryotic cell lines	<input type="checkbox"/>	<input checked="" type="checkbox"/> Flow cytometry
<input checked="" type="checkbox"/>	<input type="checkbox"/> Palaeontology	<input checked="" type="checkbox"/>	<input type="checkbox"/> MRI-based neuroimaging
<input checked="" type="checkbox"/>	<input type="checkbox"/> Animals and other organisms		
<input checked="" type="checkbox"/>	<input type="checkbox"/> Human research participants		
<input checked="" type="checkbox"/>	<input type="checkbox"/> Clinical data		

Flow Cytometry

Plots

Confirm that:

- ☒ The axis labels state the marker and fluorochrome used (e.g. CD4-FITC).
- ☒ The axis scales are clearly visible. Include numbers along axes only for bottom left plot of group (a 'group' is an analysis of identical markers).
- ☐ All plots are contour plots with outliers or pseudocolor plots.
- ☐ A numerical value for number of cells or percentage (with statistics) is provided.

Methodology

Sample preparation	Yeast cells were fixed in 70% ethanol, washed twice with 50 mM sodium citrate, treated with RNaseA (0.25 mg/ml) and ProteinaseK (1 mg/ml), and then stained with Sytox Green for at least one hour.
Instrument	MoFlo XPD cell sorter (Beckman-Coulter), BD LSRFortessa X-20 flow cytometer (Becton Dickinson)
Software	FlowJo v10.5.3
Cell population abundance	At least 5 million cells were sorted for samples intended for DNA extraction; 300000 cells were sorted as a control sample. The purity of sorted cells was assessed by restaining an aliquot of cells with Sytox green, followed by flow cytometry (see Figure 2c)
Gating strategy	Cells were initially gated using the FL1-Height and the FL1-Area signal to select single cells. A histogram of the FL1-Area was used to set the gate for replicating and non-replicating cells (see Figure 2 a)

☒ Tick this box to confirm that a figure exemplifying the gating strategy is provided in the Supplementary Information.